ACP-25

Schlafapnoe-Atemtherapiegerät - Bedienungsanleitung
Sleep apnea respiratory therapy device - Instruction Manual
Appareil de traitement de l'apnée du sommeil - Mode d'emploi
Dispositivo per la terapia respiratoria dell'apnea del sonno - Istruzioni per l'uso

www.pulox.de

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Dear customer,

we are pleased that you have chosen a product from our range. Please read these operating instructions carefully and keep them accessible for other users. Observe the application instructions and notes.

With kind regards

Your Novidion Team

ACP-25 (BR-2XXX Series)

Sleep Apnoea Breathing Therapy Equipment

Models: BR-20C, BR-20A, BR-20S, BR-20ST, BR-25C, BR-25A, BR-25S, BR-25ST.

Instruction for use

- Please read the instruction carefully before using the product.
- Please put this instruction book in a place where you can read it quickly.

1. Summary

This manual provides guidance for the installation, use and daily maintenance of BR-2XXX series Sleep Apnoea Breathing Therapy Equipment. These specifications apply to the following models: BR-20C, BR-20A, BR-20S, BR-20ST, BR-25C, BR-25A, BR-25S and BR-25ST.

Patients should read the instruction carefully before using the treatment equipment, and strictly follow the instructions to ensure safety and effectiveness.

The terms "warning" and "caution" in this specification refer to "safety" and "effectiveness" in the use of the equipment. The specific meanings are as follows:

- **Warning:** Failure to follow the warning may cause injury to the operator or patient.
- **Caution:** Failure to comply may cause equipment damage or malfunction, affecting the effectiveness of the equipment or ease of use.



Copyright of the user's manual is owned by Shenzhen Bi-rich Medical Devices Co., Ltd. Copyright is protected by copyright laws of various countries and international copyright conventions. For the above copyright content, the company reserves the right to pursue legal liability for any use other than reasonable use without the company's written permission.

1.1 Warning!

- This machine must be used under the guidance of a doctor!
- This device cannot be used as any form of life support system
- Guidelines for use in this manual are not a substitute for medical solutions.
- This equipment is prohibited from being used in nuclear magnetic resonance environment.
- It is forbidden to use this equipment in strong interference electromagnetic environment.

- If there is any abnormal treatment for unknown reasons, please immediately disconnect the power supply and stop using the device. Please contact the after-sales service personnel.
- When using, pay attention to the placement of the power cord and breathing tube to avoid hanging due to cables and hoses, especially due to excessive length.
- When changing the water tank, please avoid touching the heating plate at the bottom of the humidifier, or you may burn your skin.
- The device should be used under the guidance of doctors (such as mode selection, treatment parameters), and keep in touch with doctors during the use process.
- Please keep the unit, power cords, etc. out of reach of children.
- This device should not be used for people who are not self-conscious or who cannot remove the breathing mask by themselves.
- Patients must follow the respiratory parameters defined by the doctor to set up the equipment.
- Discontinue use if the output gas temperature is found to be thermal.
- Do not to use near or in the environment of oxygen-enriched or anesthetic gases.
- This product is not intended for use by children.
- When product needs to be used, cosmetics or skin care products should not be used to avoid chemical reactions.
- Please take good care of your device. Don't handle it rudely or place it carelessly; don't intentionally damage it or its accessories.

1.2 Symbols Description

Symbols	Description	Symbols	Description
	Direct current		Alternating current
	Caution! Please read the instructions carefully before use.		Type of protection against electric shock: Class II equipment
	Non-ionizing electromagnetic radiation		Refer to instruction manual
	Caution, hot! The maximum surface temperature is up to 75°C.		Classification of application components: TYPE BF APPLIED PART
	Serial number		Manufacturer
	Batch code		Date of manufacture
	Medical device		Unique device identifier
	Made in China		Model
	CE marking, Issued by SGS Belgium NV (NB No.: 1639)		Does not contain natural rubber latex
	Signs for the disposal of waste electrical and electronic equipment (Please comply with local laws and regulations)		The first number 2: Protection against solid foreign objects of diameter 12.5mm and greater; The second number 2: Protection against vertically falling water drops when enclosure tilted up to 150. Vertically falling drops shall have no harmful effect when the enclosure is tilted at any angle up to 150, on either side of the vertical.
	Authorized representative of European Community		Authorized representative of Switzerland

1.3 Intended use

The BR-2XXX series sleep apnea respiratory therapy device is intended for the treatment of obstructive sleep apnea (OSA) and hypopnea.

This device is designed to relieve the symptoms of patients suffering from obstructive sleep apnea by applying therapeutic breathing pressure to the patient's airway.

This device cannot be used for life support. It cannot be used to treat central sleep apnea. Patients are only allowed to set up the device under the guidance of a professional doctor and then use it. Professional doctors should ensure that patients understand the various functions and operating methods of the equipment.

1.4 Patient population

a) Age: adult

b) Patient: obstructive sleep apnoea (OSA) and hypopnoea.

1.5 Intended users

- a) Legally certified medical experts, such as doctor, nurse and therapist, or healthcare personnel or patient under the guidance of qualified medical experts
- b) The user should also be capable of understanding general operation of the device and the content of instruction manual.

1.6 Contraindications

Contraindications include absolute and relative contraindications:

(1) Absolute contraindications:

- Pneumothorax and mediastinal emphysema;
- Cerebrospinal fluid leakage;
- Craniocerebral trauma or intracranial or extracranial trauma or pneumocephalus;
- Shock caused by various causes and not corrected;
- Epistaxis active stage;
- Upper gastrointestinal bleeding not effectively controlled;
- Coma or with disturbance of consciousness cannot cooperate with or receive treatment;
- Vocal cord polyps;
- Severe bullous lung disease;
- Pathologically low blood pressure;
- Dehydration
- Children and pregnant or nursing women.

(2) Relative contraindications:

- Severe coronary heart disease with left heart failure;
- Acute otitis media; respiratory secretions and cough weakness;
- Weak spontaneous breathing (except T mode);
- Tubal intubation (through nose or mouth) and tracheotomy;
- Severe nasal obstruction caused by various reasons;
- Bullae;

Warning!

When patients suffer from severe chest pain, headache and dyspnea, they should contact the attending physician or salesperson as soon as possible.

The unit must be used with the supplied breathing tube and mask. To ensure the safety use, please make sure to use these breathing tube and Mask. Other breathing tube and Mask are not verified with this device.

1.7 Side effects

Patients should report unusual chest pain, severe headache, or increased breathlessness to their prescribing physician. An acute upper respiratory tract infection may require temporary discontinuation of treatment.

The following side effects may arise during the course of therapy with the device:


- Drying of the nose, mouth, or throat
- Nosebleed, stuffy nose
- Bloating
- Ear or sinus discomfort
- Eye irritation
- Skin rashes
- Headache, sore throat, sinus pain
- Reduced sense of smell
- Leakage (i.e. condensation of water into the PAP circuit, the patient's face and nose or mouth)
- Leaky mask

1.8 Working Principle

In patients with sleep apnea-hypopnea syndrome and severe snoring, temporary upper airway closure or partial closure may occur during sleep. In patients with respiratory insufficiency, pulmonary ventilation and/or ventilation dysfunction may occur, both of which may cause normal breathing to be interrupted or weakened, hypoxia, with or without carbon dioxide. The occurrence of retention symptoms is possible.

The treatment device uses a dedicated air compressor to provide two continuous positive pressures by utilizing the filtered air around it. When the patient breathes out, the therapeutic machine outputs low-pressure airflow; when the patient breathes in, the therapeutic machine outputs high-pressure airflow. The device transmits these two positive pressures to patients through an air tube device. The expiratory pressure keeps certain intra-airway pressure during expiration, keeps the upper airway of patients with sleep apnea-hypopnea syndrome and severe snoring open, and effectively alleviates carbon dioxide retention in patients with respiratory insufficiency. The difference of pressure increases the ventilation volume of patients with respiratory insufficiency, which ultimately improves the respiratory function of patients. The

treatment pressure required by each patient is determined by the doctor. The schematic diagram of the pneumatic system of the therapeutic machine is as follows:

- A. Air filter cotton F. Air-bleed hole
 B. Fan G. Flow sensor
 C. Pipeline H. Pressure sensor
 D. Patient I. Warm humidifier
 E. Mask

Fan: The pneumatic circuit drive device provides therapeutic pressure to the patient during operation of the device.

Pressure sensor: A pressure measuring device that can feedback real-time pressure values.

Flow sensor: A flow measurement device that can feed back real-time flow values.

1.9 Working mode

This series of equipment has the following treatment modes:

- **CPAP mode:** continuous positive pressure single horizontal output mode, continuous positive pressure ventilation according to the set Inspiratory pressure.
- **Auto CPAP mode:** According to the number of apnea or hypopnea events in a certain period of time, the pressure is intelligently adjusted between the preset treatment pressure and the maximum pressure; the pressure adjustment step range is 0.1-3 hPa; when the current treatment pressure is greater than 10 hPa, the pressure of automatic adjustment is 0.1-1.5 hPa; and the current treatment pressure is 0.1-1.5 hPa. When the treatment pressure is less than 10 hPa, the pressure of automatic adjustment is 1.6-3 hPa.
- **S mode:** The device is triggered by patient's spontaneous breathing to perform respiratory pressure conversion; if the ventilation function of the device is turned on and the patient has no spontaneous breathing within 6 seconds, the frequency of respiratory pressure conversion is 10 times per minute.
- **Auto S mode:** According to the number of apnea or hypopnea events occurring in a certain period of time, the inspiratory pressure is intelligently adjusted between the preset inspiratory pressure and the maximum inspiratory pressure. The expiratory pressure and the inspiratory pressure are adjusted synchronously, and the pressure difference remains unchanged. The range of pressure adjustment step is 0.1-3 hPa, the pressure of automatic adjustment is 0.1-1.5 hPa when the current treatment pressure is greater than 10 hPa, and the pressure of automatic adjustment is 1.6-3 hPa when the current treatment pressure is less than 10 hPa.
- **T mode:** The device can manually set the breathing frequency and inspiration time, and according to the set frequency and time, the two-level mode conversion.
- **S/T mode:** When the spontaneous breathing frequency is lower than the set breathing frequency, the device automatically converts from S mode to T mode, or from T mode to S mode when the spontaneous breathing frequency is higher than the set breathing frequency.


2. Equipment Introduction

2.1 Product Description of Different Types


Product Model	Product Description				Remarks
	Working mode	Max. working pressure	Pressure range	Display	
BR-20C	CPAP	20hPa	4~20hPa	3,5 inch	Single level
BR-20A	CPAP, Auto CPAP				
BR-20S	CPAP, Auto CPAP, S, Auto S				
BR-20ST	CPAP, Auto CPAP, S, Auto S, T, S/T				
BR-25C	CPAP	25hPa	4~25hPa		Single level
BR-25A	CPAP, Auto CPAP				

BR-25S	CPAP, Auto CPAP, S, Auto S				Double level
BR-25ST	CPAP, Auto CPAP, S, Auto S, T, S/T				


2.2 Product packing list

Name	Number	Picture	Remarks
1. Host-/Main unit	1 pc		
2. Adapter	1 pc		
3. Power cord	1 pc		
4. Heated humidifier	1 pc		
5. Instruction manual	1 pc		
6. MicroSD	1 pc		
7. Air filter cotton	2 pcs		
8. Transport bag	1 pc		
9. Product Qualification certificate	1 pc		
10. Error-Sheet	1 pc		


2.3 Product structure

Name	Functional description
A. Air Outlet	Connect directly to the tube
B. Guide column	When the main engine and the heater and humidifier rotate and dock, they assume an important role.
C. Host and humidifier lock	Assumes the role of the locking mechanism when the main motor and heater unit rotate and dock with the humidifier.
D. Humidifier docking ring	When the air outlet of the main motor is connected to the heating unit and the humidifier, it assumes the role of rotating and blocking the centre while preventing air from leaking out.
E. Contact area for contact pins	Connection for heating and humidifier
1-6: Locks	For the connection of the host and the heated humidifier
7-10: Contact pins	Connect to the Contact area of the heated humidifier
F. Display screen	Display operation menu, prompt and monitor data and other information


Name	Functional description
A. Heater and humidifier lock button	After pushing the button to the left, you can revolve and apart the host and the humidifier.
B. Buzzer hole	Makes the buzzer sound.
C. TF card socket	Please insert 8GB TF card. It is used for data storage.
D. USB socket	Used for the programming of the main MCU
E. Power supply socket	DC 005 socket, 24V DC power plug.

Name	Functional description
A. Water tank window	Through this window, you can see the water level of the water tank. When the water is heated, the LED lights in the water tank window light up.

2.4 Product Description

Name	Functional description
A. Power Standby 	Shortly press this button to start ventilation. Long press can stop the function. In standby mode, press the [power standby] and [Select] keys for more than 3 seconds to enter the setup interface of the treatment equipment. In the ventilated state, press the [power standby] and [Select] buttons for more than 3 seconds to enter the parameter setting interface of the current mode of operation.
B. Selection/confirmation key 	Selection and confirmation of menu items
C. Return key 	Return to the menu above or to the main interface.
D. "▲" "▼" key	It is used for cursor movement and setting parameters or item changes.

2.5 Initial use of equipment

2.5.1 Placement

Horizontal placement of equipment on a safe and reliable platform (correct placement and positioning of equipment is critical to the continuous operation of equipment).


Warning!

- When the equipment is used, it must be far away from water sources.
- This equipment should not be used in high humidity environment and liquid is not allowed to enter the equipment.
- Once the equipment falls off, the case is damaged or the water enters, please stop using immediately and contact the customer service personnel.
- This equipment must be used in the environment below 35°C and above 0°C. If the temperature exceeds this limit, it may cause irritation to the human airway.
- If the equipment is stored in an overheated or supercooled environment before use, it should be placed 2 hours in advance. In normal temperature environment, it can be reused only if the indoor temperature in which it is used is the same.
- This appliance should not be placed near heating or cooling appliances (such as radiators or air conditioners).
- It is not allowed to block the ventilator of the equipment with sheets, quilts, towels and pillows, so as to avoid blocking the ventilator.
- It is strictly forbidden to use this equipment in flammable and explosive environment (such as anesthetic gas environment and high oxygen concentration places), otherwise it may cause fire and burns.
- This device is only for personal use. To avoid the risk of cross infection, please do not share your accessories with others.
- Smoking should be avoided in the environment where the equipment is used. Tobacco smoke accumulates in the equipment, resulting in abnormal operation of the device.
- The known conditions that will adversely affect the performance of the equipment are: exposure to electric cautery, electrical surgery, defibrillation, X-ray (gamma ray), infrared radiation, transient magnetic field including magnetic resonance imaging (MRI), radio frequency interference.

2.5.2 Installing

(1) Installation of air filter cotton:

Install the air filter cotton and the air intake cover on the main unit as shown in the figure below.


- A. Air filter cotton
B. Intake cover

Warning!

Before use, it is necessary to ensure that air filter cotton is installed in the equipment.

(2) Connecting power lines

- a.) Connect the plug of the power adapter to the device;
- b.) Connect the power cord to the power adapter;
- c.) Insert the power cord plug into the socket.


A. Power cord

Note: When the power supply is interrupted and restored, the equipment will automatically return to the working state before the power supply is interrupted.

⚠ Caution!


- Always check the power cord for damage, as damage to the power cord will cause the unit to not function properly.
- It is necessary to check whether the power cord is damaged or not. Damage to the power cord will cause the unit to stop working properly.

(3) Connect the tube to the equipment and insert the tube into the outlet as shown below.

A. Air tube

(4) Connect the end of the tube to the mask interface as shown below.

A. Mask
B. Air tube

(5) Insert the microSD card: Insert the microSD card into the microSD socket

2.5.3 Start treatment

In standby mode, press the [power] button to start treatment.

⚠ Caution!

- In the ventilation mode, long press the [power] button to stop ventilation.
- In order to prevent the disconnection between the breathing tube and main unit or between the breathing tube and mask during operation, only pipes conforming to ISO 5367 or ISO 80601-2-74 should be used.
- The maximum pressure of the equipment in use is:
For BR-20C/ BR-20A/ BR-20S/ BR-20ST: **20hPa**
For BR-25C/ BR-25A/ BR-25S/ BR-25ST: **25hPa**

Warning!

- When the patient has severe chest pain, headache, chest distress and dyspnea, please stop the treatment and contact the attending physician at the first time.
- Patients must follow the respiratory parameters defined by the doctor to set up the equipment.
- Covering the tube with blankets or heating equipment with insulating boxes may affect the treatment results or cause injury to patients.
- Not using masks or accessories that minimize the rebreathing of carbon dioxide or allow spontaneous breathing can cause asphyxiation.

- The accuracy of leak alarm will be affected if the equipment is used when atmospheric pressure exceeds the range (see Part 9.1, Operating Environment).

2.6 Daily use

2.6.1 Connect tube, wear mask and headband

Apply the mask and the headband. Please follow the instructions below during this process.

⚠ Caution!

- The packaging integrity of disposable breathing tube should be checked before use. If the package is broken, please do not use it and replace it with a new one.
- Check whether the mask is clean before use, if not, clean the mask according to the instruction for use of mask provided by the mask manufacturer.
- Do not stack excessively long tubes at the head of the bed, possibly entangling the head or neck while the patient is asleep. The tube should not be covered or close to the heating source, which may cause dangerous deformation of the tube.

2.6.2 Adjustment of Tube

The patient can adjust the tube after lying in bed, so that the tube can also move freely when the patient turns over in a sleeping state. Adjust the mask and headband airtight to make the patient feel comfortable.

2.6.3 Initiation ventilation

Shortly press the button [power standby] after booting, the device starts to work, and the display screen will display the relevant treatment information.

2.6.4 Use Boost Delay Function

After selecting the boost delay function, the output pressure of the equipment will gradually rise from the initial pressure and reach the set treatment pressure within the preset delay time. When the pressure rises gradually, the patient can fall asleep easily. The display will show the remaining time of the delay boost.

2.6.5 Shut down the equipment

Remove the mask and headband, press the [power] button for 3 seconds, and the device will stop working.

2.7 Use of Heating and Humidifier Device

Temperature humidifier and ventilator are used at the same time, which will not change the expected use of the ventilator.


Nasal discomfort can easily occur if the ventilator is used for a long period of time in an area with dry air or at night. The device can increase the humidity of the inhaled air, protect the nasal mucosa from drying out and improve comfort.

WARNING:

- To preserve your health, you must observe the requirements of 5.2.2 for disinfection before each use of the humidifier!

2.7.1 Water is injected into the heater and humidifier tank

(1) Water is injected directly from the top of the tank to the maximum water level calibration line, with a water addition of about 270mL, which should not exceed the maximum water level.


A. Kettle
B. Water tank

Precaution!


- Be sure to empty the remaining water before adding water.
- Pure water or cold boiled water should be added to the water tank, and no other substances should be added. Otherwise, it may not be able to be used normally, affecting the quality of treatment or injuring patients.

2.7.2 Pour out the water in the tank

(1) Separate the main engine from the heater and humidifier first. As shown in the figure below, pull the push button to the right and rotate the heater and humidifier counterclockwise.


(2) Then fix the heater and humidifier with the right hand, press the lock button of the water tank with the left hand, gently pull the water tank out, and separate it from the heater and humidifier.

1:


- A. Heated humidifier
- B. Water tank lock button

2:

- A. Heated humidifier
- B. Water tank

(3) Take out the water tank, open the lock of the water tank and open the water tank.

- A. Water tank buckle
- B. Water tank

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⚠ Caution!

When not using the heater and humidifier, empty the water inside the tank and dry it.

2.7.3 Installation of Heating and Humidifier Device

(1) Select "Settings" and then "Heating and Humidifier Device".

(2) When "off" is selected, the heater and humidifier will not work.

(3) Adjust gear according to "▲" and "▼". there are five gears, from 20% to 100%, the working temperature of the heater and humidifier increases gradually.

⚠ Caution!

- Generally, when the temperature in the water tank is low, the humidity of the air in the mask is also low.
- When the indoor temperature is too low and the temperature of the tube is too low, the content of the tube is prone to condensation.
- If the patient's nasal cavity feels dry, indicating that the humidity is too low, the warming and humidifying gear should be raised.
- If significant air volume change is found, stop using and change the device


Warning!

- Please use this device under the specified environmental conditions, beyond the operating conditions may cause device malfunction or injury to the patient.
- When using the heater and humidifier, observe the water level in the water tank to ensure that there is enough water in the water tank (the maximum water level must not be exceeded by refilling water) and do not heat the heater and humidifier dry.
- When the heater and humidifier work, do not touch the heating plate to avoid scalding.
- It is required to pour out the excess water after use and dry the tank.

2.8 Carry treatment equipment

(1) Take off the air hoses and empty the water in the water tank according to the method described in 2.7.2.

(2) Pack all pieces of treatment equipment in the transport bag.

- | | |
|------------------------|----------------------|
| A. Host | D. Heated humidifier |
| B. Accessories package | A. E. Transport bag |
| C. Air tube | |

Warning!

Before carrying the equipment, the water tank must be emptied to prevent the remaining water from entering the equipment.

3. Patient menu parameter setting

3.1 Parameter Setting Steps

(1). Connect the power supply and the power adapter. The display screen displays the main page as follows:

(2). Press and hold the [Power Standby] and [Select] keys at the same time for 3 seconds, then the parameter setting interface can be accessed. The cursor will switch between the icons, and the screen will display different interfaces. As shown in the following figure:

(3). Press [▲] [▼], when the cursor stays at the option that needs to be adjusted, press the button to enter the adjustment state. At this time, the color of the background of the option is green, and the adjustable font is red. As shown in the following figure:

(4). Press [▲] [▼] to adjust the working parameters, when the setting value is selected, press the [Select] key, then the setting value is determined, and the background color of the option becomes blue. As shown in the following figure:

(5). When the settings are completed, press the [Return] button, and the parameter settings interface will fall back to the original interface.

3.2 Operating Interface Description

The device provides four operation interfaces, namely "operation mode setting interface", "user information display interface", "user setting interface" and "factory setting interface".

- "Operating mode setting interface", i.e. working mode parameter setting interface: under this interface, there are six operating modes: "CPAP", "Auto CPAP", "S", "Auto S", "T" and "S/T". Different parameters must be set for each operating mode. See 3.3 for details.
- "User Information Display Interface": In this interface, besides setting up "statistical cycle", other information is provided to users for device operation.
- "User Setting Interface": in this interface, you can set time, time format, unit, heating humidifier and so on.
- "Factory Setup Interface": It's not for users, there's no explanation.

3.3 Detailed description of setting options and menus

Make a setting	Detailed description	Applicable model
Pattern	There are six operation modes: CPAP mode, Auto CPAP mode, S mode, automatic S mode, T mode and S/T mode.	/
CPAP	"Continuous positive airway pressure" (CPAP) is based on therapeutic pressure and continuous positive airway pressure ventilation.	All models
Auto CPAP	It can automatically detect the patient's apnea and airflow reduction, and then according to the pressure of upper respiratory tract resistance, sleep phase,	Except BR-20C BR-25C

	position of different pressure. With the minimum output pressure to achieve the best treatment effect and better comfort.	
S	It is also known as "independent trigger mode" or "synchronization mode". Triggered breathing is a person's own breathing. The working frequency of this mode is completely controlled by the patient's own breathing. It is mainly suitable for patients with good breathing.	BR-20S BR-20ST BR-25S BR-25ST
Auto S	Based on the S model, the inspiratory pressure can only be adjusted between the preset inspiratory pressure and the maximum inspiratory pressure according to the number of apnea or hypopnea events in a certain period of time. The expiratory pressure and inspiratory pressure are synchronously adjusted with the pressure difference unchanged.	BR-20S BR-20ST BR-25S BR-25ST
T	The device will be controlled by two-stage conversion according to inspiration time and breathing frequency.	BR-20ST BR-25ST
S/T	When the apnea time reaches 10 seconds, the device will automatically switch from S mode to T mode. Once automatic inhalation is detected, the device will switch back to S mode.	
CURA pressure	"Curative Pressure Applicable" for CPAP, Auto CPAP mode. The therapy equipment continues to provide the patient with the pressure of the ventilation.	All models
INIT Pressure	Initial pressure: The pressure will be used as the curative pressure when the therapy equipment starts. Then the curative pressure was gradually increased to the set of curative pressure according to the remaining ramp time.	All models
Ramp Time	0~60 minutes, 5 minutes each step. During the ramp time, the curative pressure will gradually increase from "initial pressure" to the set "curative pressure".	All models
Manual P Adj	Manual pressure adjustment 0~2hPa, 0.5hPa each step. When the therapy equipment operates in this mode, the [▲] and [▼] button can be used to regulate the curative pressure in real time. The range of pressure regulation is the curative pressure. Press ▲/▼ to manual adjust the value.	All models
MINI Pressure	The "minimum curative pressure" applicable to Auto CPAP mode. In Auto CPAP mode, the therapy equipment will operate in this curative pressure after the ramp time.	Except BR-20C BR-25C Except BR-20C BR-25C
MAX Pressure	The "maximum curative pressure" applicable to Auto CPAP mode. In Auto CPAP mode, the therapy equipment will operate in this curative pressure after the ramp time. Then the therapy equipment will adjust the curative pressure intelligently according to the number of apnea or hypoventilation events occur within a certain period of time.	
INSP Pressure	Inspiratory pressure applicable to S mode and Auto S mode. The pressure set for the inspiratory process by the therapy equipment.	BR-20S BR-20ST
EXPI Pressure	Expiratory pressure applicable to S mode and Auto S mode. The pressure set for the expiratory process by the therapy equipment.	BR-25S BR-25ST
Inhalation Time	Inspiratory time applicable to T mode and S/T mode. Range: 0.5~3 seconds	BR-20ST
RESP Rate	Respiratory rate applicable to T mode and S/T mode. Range: 6~40	BR-25ST
Maximum Time	The maximum inspiratory time Range: 0.3~3 seconds applicable to S mode. The real inspiratory time should be no more than the maximum inspiratory time.	BR-20S BR-20ST BR-25S BR-25ST
Minimum Ti	The minimum inspiratory time Range: 0.3~3 seconds applicable to S mode. The real inspiratory time should be no less than the minimum inspiratory time.	
pressure Rising Time	"Pressure rising Time": 1-4 gear adjustable, applicable to S mode and Auto S mode. The greater the value, the longer the pressure rising time.	
Adjusting pressure sensitivity	Pressure adjustment sensitivity: three options (low, medium, high). Only applicable to Auto CPAP mode. When the option is high, the pressure adjustment sensitivity is the highest one.	Except BR-20C BR-25C
Inspiratory sensitivity	Inspiratory sensitivity Range: 1-8. The larger the setting value, the higher the trigger sensitivity and the smaller the trigger flow.	BR-20S BR-20ST
EXPI sensitivity	Expiratory sensitivity Range: 1-8. The larger the setting value, the higher the trigger sensitivity and the larger the trigger flow.	BR-25S BR-25ST
Backup VENT	Backup Ventilation: Only applicable to S mode. When the backup ventilation function is turned on, the therapy equipment will start 10 times per minute backup ventilation when the patient does not breathe spontaneously in 6 seconds. When the backup ventilation function is closed, there is no function.	BR-20S BR-20ST BR-25S BR-25ST

MAX INSP	The maximum inspiratory pressure: only applicable to Auto S mode. The curative pressure will be automatically adjusted between the inspiratory pressure and the maximum inspiratory pressure.	BR-20S BR-20ST BR-25S BR-25ST
Heater	Heated humidifier: 20%-100% adjustable. The greater the gear, the higher the temperature and humidity.	All models
Expiratory pressure drop	When the patient exhales, the curative pressure is automatically decreased. "ON" indicates that the function is enabled. "Off" indicates that the function is disabled.	All models
Time Setting	Through [▲], [▼] and [Select] button, you can set up the time in detail.	All models
Leakage Prompt	Once the leakage volume exceeds 85 LPM/min for more than 36 seconds, a leakage prompt will appear on the title bar.	All models
Auto Power On	When the auto power on function is on, once the airway pressure variation in one second is more than 1 hPa, the therapy equipment will start automatically with the last operating parameter. When the auto power on function is off, there is no such function	All models
Auto Power Off	When the automatic power-off function is activated, the treatment device will automatically stop ventilation if the air leakage exceeds 65 L/min for more than 3 minutes. The auto power off function is not available when it is turned off.	All models

4. Buzzing Warning

4.1 Beep warning level description

Grade	standard	description
High priority	***	Operators must respond immediately.
Medium priority	**	Operators must respond.
Low priority	*	The state of the equipment has changed, which requires the operator's caution.

4.2 Buzzing Warning

Alarm level	Sound Description
High priority	dididi didi dididi didi
Medium priority	dididi dididi
Low priority	di di

4.3 Buzzing Warning Types

Type	Warning Priority	Description
Power Error!	high	Power Error: If the power supply is not normal, the microcontroller cannot work properly, and may not be able to send an effective alarm.
High temperature of motor!	high	Motor temperature: When the switch MOSFET is pierced, the temperature of the motor may be too high.
Leak!	middle	Leakage: Leakage exceeds 65 L/min, lasting more than 36 seconds.
Low pressure!	middle	Low pressure: The treatment pressure is less than 3 hPa and lasts for more than 30 seconds.
SD card error!	middle	MicroSD card error: If the SD card cannot be written, an alarm will be generated.

5. Cleaning and maintenance


Warning!

- Regular cleaning of equipment accessories can effectively prevent respiratory infections.
- Please always purchase breathing accessories that meet the requirements and are marked with CE mark. The use of unqualified accessories may reduce the safety of medical devices.
- Be sure to use detergents that are harmless to the human body and do not cause allergic reactions.
- Be sure to remove the power plug before you start cleaning, otherwise it will cause an electric shock.
- The equipment should not be immersed in any liquid, and this should also be avoided during cleaning.
- Do not touch the heater before it cools, otherwise it may damage the skin.

5.1 Daily cleaning

Cleaning water tank

- (1) Open the water tank according to the method of 2.7.2, immerse the soft cloth in a little detergent, wipe the top and bottom of the water tank, then rinse it with tap water, then wipe it with soft cloth, and place it in the dry and ventilated place to dry naturally.

A. Upper cover

B. Lower tank

(2) After the water tank is dried, cover the water tank and put the water tank back into the heater and humidifier.

5.3 Weekly cleaning

Recommended cleaning agent: 70% alcohol.

5.2.1 Cleaning

- Use gauze to dip in soapy water (no more than 40°C). Wipe the surface of the main machine and the heater and humidifier.
- Water tank clean procedure: Refer to the clause 5.1.

Warning! The equipment can only be used after the case is dried.

5.2.2 Disinfection

If you follow the above cleaning requirements throughout their service life and clean well, it is generally not necessary to disinfect your equipment. When the equipment is contaminated or used in clinic, it must be disinfected according to the following methods.

3. Disinfection of surface of the main machine and the heater and humidifier:

(1) Use neutral disinfectant (70% alcohol).

4. Disinfection of water tank of heating and humidifying apparatus:

Cleaning before disinfection is carried out in accordance with Part 5.1 "Cleaning water tank".

The disinfection procedure as below:

Step 1: Spray the inside and outside of the water tank with 70% alcohol, then wipe with a clean cotton towel;

Step 2: Immersed in tap water for 30 minutes at a temperature of 75°C +2°C;

Step 3: After taking out, drain off the water and wipe dry with a dry cotton towel.

5.3 Replacement of air filter cotton

(1) Open the intake cap according to 2.5.2 and remove the air filter cotton.

(2) Put the new air filter cotton into the filter, and then cover the air cap.

⚠ Caution!

- To avoid damage, spare air filter cotton should not be placed under direct sunlight, humid environment and freezing temperature.
- Air filter cotton should be replaced at least every 6 months. If it is obviously dirty, damaged or cracking, it should be replaced in time.

6. Maintenance

⚠ Caution!

- Maintenance and maintenance of treatment equipment should be carried out by users.
- Under no circumstances may the user open the housing of the treatment unit. Repair and maintenance of the indoor unit may only be carried out by the supplier.
- If you have abnormal phenomena in the use process, such as unstable pressure, you need to contact after-sales service personnel for testing and maintenance.

6.1 Warranty Conditions

a.) Free maintenance can be carried out during the warranty period due to faults caused by equipment materials or manufacturing problems.

b.) The warranty period shall be calculated on the basis of the information on the repair card or from the date of purchase invoice of the equipment.

c.) After-sales service hotline: 0755-27464589

6.2 Storage

Always store your medical devices in a dry and ventilated environment. Do not store them in a dark or damp environment!

6.3 Disclaimer

If there are the following circumstances, the company will no longer bear any responsibility for the product:

- a.) Damage caused by disassembly and repair by maintenance personnel other than our company.
- b.) Damage caused by failure to comply with these instructions.
- c.) Damage to equipment caused by the use of accessories not approved by the company.

7. Waste disposal

⚠ Please ensure that the equipment and accessories other than its service life are disposed of in accordance with the relevant laws and regulations of the State. If the equipment is not disposed of in accordance with the requirements, it will cause harm to the environment. Generally, equipment with end-of-service life should be sent to the recycling mechanism with cardboard or protective plastics. The recycling mechanism should be able to dispose of plastics, glass, metal parts, printed circuit boards, wires and cables, heating boards of water tanks and motor materials.

8. Common problems in use and troubleshooting

The following are common problems and solutions in the use of some equipment. If the following solutions fail to achieve the desired results, contact your doctor or equipment supplier.

8.1 Common problems and solutions of patients:

Problem	Possible causes	Terms of settlement
Dry mouth and throat	Non-heating humidifier	Turn on the heater and humidifier
	Temperature of heater and humidifier is too low	Increase the temperature of heating and humidifying device
Red face	The mask is too tight	The mask should be as loose as possible, without air leakage.
	The size of the mask is not right	Everyone's face is different, you may consider changing the mask.
	Allergy to mask materials	Consult a doctor or equipment supplier Place a skin-friendly material between the skin and the mask.
Dyspnea	Treatment pressure or excessive expiratory pressure setting	Consult a doctor to reduce the pressure of treatment appropriately.
Difficulty in inhalation	Air leakage and rupture of tube	Contact the after-sales service personnel to replace the tube.
	Treatment pressure or inhalation pressure setting too low	Consult the doctor and adjust the treatment pressure properly.
	Inhalation sensitivity set too low	Please choose high sensitivity inhalation.
	The intake is blocked.	Remove the object blocking the intake.
Pain in the nose or ear	Inflammation of sinus or middle ear	Contact your doctor immediately
Unable to adapt to treatment stress, feeling uncomfortable	When the treatment pressure is higher than 13 hPa, patients will feel uncomfortable, but the pressure is set according to the condition, too low cannot prevent apnea.	Patients will have a process of adaptation to pressurized air for up to 4 weeks. They can try to relax themselves and learn to breathe through their nose. If the problem persists, please contact your doctor.

8.2 Common Faults of Equipment and Solutions

Problem	Possible causes	Terms of settlement
Mask pressure is different from actual pressure	Improper tube connection	Connect the tube to the correct interface.
	A leak in a mask or pressure sensor	Please contact the equipment supplier.
Equipment does not work	Automatic power off function enabled	Press the power button to restart the device.
	Incorrect insertion of power plug	Plug in the power plug correctly.
	Poor contact of power plug	Please find a professional electrician to replace the power socket.
Standby status cannot start	Procedures for treatment equipment did not start properly	Disconnect the power supply of the device, wait 30 minutes before inserting the power supply, and then restart.
Excessive noise	Improper tube connection	Connect the tube to the correct interface
Output temperature of sensory equipment is abnormal	The air path of the unit may not be smooth and the incoming air may not be sufficient.	Replace air filter cotton and clean vent
		Keep away from objects that interfere with air circulation
Low output gas pressure	Airway may be blocked	Change the air filter cotton and clean the air vent to check whether the intake passage of the equipment is unblocked.
	If the delayed boost function is used, it will take a certain time for the pressure to rise to a certain level. Normal, if the time is within	If necessary, the delay boost function can be turned off or a shorter delay boost time can be set.

the range.

9. Technical Specifications of Products

9.1 Equipment Performance Index

For the unit																																
Treatment mode	BR-20C and BR-25C	CPAP																														
	BR-20A and BR-25A	CPAP and Auto CPAP																														
	BR-20S and BR-25S	PAP, Auto CPAP, S and Auto S																														
	BR-20ST and BR-25ST	CPAP, Auto CPAP, S, Auto S, T, S/T																														
Continuous running time	≥8 hours																															
Maximum inspiration time	3s																															
Operating pressure range	For BR-20C/ BR-20A/ BR-20S/ BR-20ST: 4-20 hPa (cm H ₂ O); For BR-25C/ BR-25A/ BR-25S/ BR-25ST: 4-25 hPa (cm H ₂ O);																															
Pressure error (between the measured value and the set value)	± 0.5 hPa																															
Dynamic Pressure Accuracy	+0.5/-2.0hPa.																															
Maximum single fault steady state pressure	For BR-20C/ BR-20A/ BR-20S/ BR-20ST: 30 hPa; For BR-25C/ BR-25A/ BR-25S/ BR-25ST: 40 hPa																															
Maximum gas flow	<p>When the Operating pressure is set to the values listed in the following table, the average flow rate of the patient interface should be greater than 80% of the corresponding flow value in the following table. The following table is applicable to BR-20A, BR-20C, BR-20S, BR-20ST (tested in CPAP mode)</p> <table border="1"> <tr> <td>Set pressure (hPa)</td> <td>4</td> <td>10</td> <td>15</td> <td>20</td> </tr> <tr> <td>Pressure measured at the patient's connection (hPa)</td> <td>3</td> <td>9</td> <td>14</td> <td>19</td> </tr> <tr> <td>Average patient interface flow (L/min)</td> <td>60</td> <td>65</td> <td>65</td> <td>60</td> </tr> </table> <p>The following table is applicable to BR-25A, BR-25C, BR-25S, BR-25ST (tested in CPAP mode)</p> <table border="1"> <tr> <td>Set pressure (hPa)</td> <td>4</td> <td>11</td> <td>18</td> <td>25</td> </tr> <tr> <td>Pressure measured at the patient's connection (hPa)</td> <td>3</td> <td>10</td> <td>17</td> <td>24</td> </tr> <tr> <td>Average patient interface flow (L/min)</td> <td>60</td> <td>63</td> <td>60</td> <td>60</td> </tr> </table>		Set pressure (hPa)	4	10	15	20	Pressure measured at the patient's connection (hPa)	3	9	14	19	Average patient interface flow (L/min)	60	65	65	60	Set pressure (hPa)	4	11	18	25	Pressure measured at the patient's connection (hPa)	3	10	17	24	Average patient interface flow (L/min)	60	63	60	60
Set pressure (hPa)	4	10	15	20																												
Pressure measured at the patient's connection (hPa)	3	9	14	19																												
Average patient interface flow (L/min)	60	65	65	60																												
Set pressure (hPa)	4	11	18	25																												
Pressure measured at the patient's connection (hPa)	3	10	17	24																												
Average patient interface flow (L/min)	60	63	60	60																												
Pressure drop caused by heating of humidifier	Maximum value: 0.4 hPa at 60 L/min flow																															
Water tank capacity of humidifier	270 mL																															
Humidifying rate	12.18mg/L																															
Heating temperature for humidifier	Maximum 75°C																															
Boost time	Adjustable range: 0-60 minutes.																															
The maximum inspiratory time	Should not exceed 3 seconds to prevent hyper-ventilation of a patient.																															
Leak indicator	Exceeds 65 L/min																															
Automatic air leakage compensation	When the air leakage is less than or equal to 65L/min, the pressure error ≤±1 hpa.																															
Noise	Less than 45 dB A (according to ISO 80601-2-70)																															
Classification of application components	Type BF																															
Grade of waterproof	IP22																															
Device life	5 years																															
Size (L * W * H)	300*162*107mm																															
Gross weight	1.75KG (no water)																															
Air outlet	22 Millimeter Interface (outer diameter)																															
For the accessories																																
SD card	Type: Micro SD																															
	Capacity: 8 GB																															
	Storage parameters: Treatment of pressure, flow, insufficient ventilation per hour																															
Adaptor	Input: 100-240V~; 50HZ/60HZ; 1.6 A;																															
	Output: 24Vd.c./ 2.7A																															
	Type of protection against electric shock: Class II																															
Operating environment																																
Temperature	0°C ~ 35°C																															
Humidity	≤ 93%Rh, Non-condensation																															
Atmospheric pressure	700 hPa ~ 1060 hPa																															
Storage environment																																
Temperature	-20°C ~ 60°C																															
Humidity	≤93%Rh, non-condensation																															
Atmospheric pressure	700 hPa ~ 1060 hPa																															

Statement:

All breathing air pathways of the BR-2XXX Series Sleep Apnoea Breathing Therapy Equipment were meet the ISO 18562-2 (particulate matter) and ISO 18562-3 (VOCs).

Warning!

Do not use a warming humidifier beyond the maximum rated height or beyond the rated temperature range. The use of a

warming humidifier outside this temperature range will affect the quality of treatment or harm patients.


9.2 Inspiratory and expiratory resistance

Autonomous breathing when the equipment is in danger of normal electrification due to power failure or partial power failure:


Flow Settings (L/min)	Inspiratory and expiratory resistance (hPa)
30	≤ 1
60	≤ 3

9.3 Pressure capacity curve


1/3 Maximum pressure (unit: hPa)

2/3 Maximum pressure (unit: hPa)

3/3 Maximum pressure (unit: hPa)

Note: The above parameters are measured under ATPD environment. The manufacturer reserves the right to change the parameters.

Warning:

When the equipment runs beyond the scope specified in the above 9 chapters, it may lead to testing of the consequences of inaccurate quantity or output.

9.4 Display parameter table

Parameter	Definition
P (Pressure waveform)	Respiratory pressure waveform changes with time. Average the past five pressure samples and refresh every 200 ms
F (Flow waveform)	Respiratory flow waveform changes with time, real-time display sampling every 200 ms
Pressure	Output pressure real-time monitoring value, average the past five pressure sampling, refresh when pressure change is 0.1 or more than 100 ms.
Respiratory frequency	Respiratory frequency detection value, average of the past 5 respirations
Minute ventilation	Minute ventilation, average of the past five breaths
Inhalation time	Inhalation time, last breathing time
Air leakage	Leakage detection value, average of the past five breaths, refresh every second when changing.

10. EMC Information

The equipment needs special information on EMC and the EMC information provided in this chapter.

⚠ Caution!

- The equipment of this department meets the requirements of IEC60601-1-2 standard for electromagnetic compatibility
- Users should install and use the EMC information provided by random files.
- Portable and mobile RF communication equipment may affect the performance of the equipment, avoid strong electromagnetic interference when used, such as close to mobile phones, microwave ovens, etc.
- Detailed instructions and manufacturer's statements are attached.

Warning!

- This equipment should not be used close to or overlay with other equipment. If it must be used close to or overlay, it should be observed to verify whether it can operate normally under the configuration used.
- Class A equipment is intended to be used in industrial environments. Due to the conduction and radiation disturbances of equipment, there may be potential difficulties in ensuring EMC in other environments.
- In addition to transducers and cables sold by the manufacturer of the equipment or system as spare parts for internal components, the use of accessories, transducers and cables other than those specified may result in increased emission or reduced immunity of the equipment or system.
- The use of non-prescribed accessories, transducers or cables in conjunction with equipment and systems may result in an increase in equipment or system emission or a decrease in immunity.

Table 1: Electromagnetic Emission

Guidance and manufacture's declaration – electromagnetic emission		
This is intended for use in the electromagnetic environment specified below. The customer of the user should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	Use RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.

RF emission CISPR 11	Class B	Is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Table 2: EMI

Guidance and manufacture's declaration – electromagnetic immunity			
This is intended for use in the electromagnetic environment specified below. The customer or the user should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines; ±1 kV for input/output lines	±2kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ; ± 2 kV line(s) to earth	±1 kV differential mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T) for 0.5 cycle, 40% U_T (60% dip in U_T) for 5 cycles, 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	<5% U_T (>95% dip in U_T) for 0.5 cycle, 40% U_T (60% dip in U_T) for 5 cycles, 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50Hz/60Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE U_T is the a.c. mains voltage prior to application of the test level.			

Table 3: EMI

Guidance and manufacture's declaration – electromagnetic immunity			
This is intended for use in the electromagnetic environment specified below. The customer or the user should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	3 V _{rms}	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2\sqrt{P}$ $d = 1.2\sqrt{P}$ 80 MHz to 800 MHz $d = 2.3\sqrt{P}$ 800 MHz to 2.5 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device			
b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

Table 4: Recommended Safety Distance

Recommended separation distances between portable and mobile RF communications equipment and the device.			
This is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 KHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz

0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d inch (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.


NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Appendix maintenance menu and service menu parameter setting

1. Maintenance menu parameter setting steps

This menu is for clinicians only.

<p>Enter the main interface Connect the power cord to the power adapter and display the main interface as shown on the right.</p> <p>Applicable models: BR-20C, BR-20A, BR-20S, BR-20ST, BR-25C, BR-25A, BR-25S, BR-25ST</p>	
<p>2. Setting the initial interface by setting parameters</p> <p>After entering the main interface, hold down the [power standby] key and [Select] key at the same time for 3 seconds, to access the parameter setting interface. The cursor switches between the icons, and the screen will display different interfaces. As shown on the right:</p>	
<p>3. Enter the parameter setting interface</p> <p>When the cursor stays on the icon, press the [Select] key, to enter the parameter setting interface. At this time, the first option of parameter setting, the background color becomes green, as shown on the right:</p>	
<p>Select the options to set parameters</p> <p>Press [▲] [▼], when the cursor stays at the option that needs to be adjusted, press the button to enter the adjustment state. At this time, the color of the background of the option is green, and the adjustable font is red.</p>	
<p>5. Changing option parameters</p> <p>Press [▲] [▼] to change the working parameters.</p>	
<p>6. Determine option parameters</p> <p>When the setting value is selected, press the [Select] key, then the setting value is determined. The background color of the option becomes blue.</p>	

7. Return parameter settings initial interface

Press the [Return] key, and the parameter settings interface falls back to the original interface.

Main page 5/8/2019 9:10PM

Pressure:

0.0 hPa

Ramp:

0 min

2. Maintenance menu setting options and their description

This menu is for clinicians only.

Function	Does the following pattern have this function?						Set range	Description
	CPAP	Auto CPAP	S	Auto S	T	S/T		
Therapeutic stress	Yes	No	No	No	No	No	4~20hPa	This is set to the fixed output pressure after the delay boost. Every press on [▲] or [▼] increases or decreases 0.5 hPa. The default value is 6 hPa
Minimum therapeutic pressure	No	Yes	No	No	No	No	4~20hPa	This is set to the fixed output pressure after the delay boost. Every press on [▲] or [▼] increases or decreases 0.5 hPa. The default value is: 4 hPa
Manual Voltage Regulation	Yes	No	No	No	No	No	0~2 hPa	This setting allows the patient to manually adjust the pressure during treatment by pressing [▲] or [▼] keys once, for increasing or decreasing the value by 0.5 hPa. The default value is: 0 hPa
Adjusting pressure sensitivity	No	Yes	No	No	No	No	low/middle/high	Press the [▲] or [▼] button to select the mode. The default value is: middle
Boost delay	Yes	Yes	Yes	Yes	Yes	Yes	0~60 Minute	The maximum boost time is 0-60 minutes. Every time you press [▲] or [▼] the value increases or decreases by 5 minutes. The default value is 0 minutes.
Initial pressure	Yes	Yes	Yes	Yes	Yes	Yes	4~20hPa 4~25hPa	Every time you press [▲] or [▼], the value increases or decreases by 0.5 hPa. The default value is 4 hPa
Maximum treatment pressure	No	Yes	No	No	No	No	4~20hPa	Every time you press [▲] or [▼], the value increases or decreases by 0.5 hPa. The default value is: 15 hPa
Maximum Inspiratory pressure	No	No	No	Yes	No	No	4~25hPa	Every time you press [▲] or [▼], the value increases or decreases by 0.5 hPa. The default value is: 15 hPa
Expiratory pressure	No	No	Yes	Yes	Yes	Yes	4~20hPa 4~25hPa	Every time you press [▲] or [▼], the value increases or decreases by 0.5 hPa. The default value is: 4 hPa
Suction pressure	No	No	Yes	Yes	Yes	Yes	4~20hPa 4~25 hPa	Every time you press [▲] or [▼], the value increases or decreases by 0.5 hPa. The default value is: 12 hPa
Ti max	No	No	No	Yes	Yes	No	0.1~3.0 S	Every time you press [▲] or [▼], the value increases or decreases by 0.1 S. The default value is: 3.0S
Ti min	No	No	No	Yes	Yes	No	0.1~3.0 S	Every time you press [▲] or [▼], the value increases or decreases by 0.1 S. The default value is: 0.1S
Pressure rise time	No	No	Yes	Yes	Yes	Yes	1~4	Every time you press [▲] or [▼] the value increases or decreases by 1. The default value is: 2
Inspiratory triggering sensitivity	No	No	Yes	Yes	Yes	No	1~8	Every time you press [▲] or [▼], the value increases or decreases by 1. The default value is: 6
Expiratory triggering sensitivity	No	No	Yes	Yes	Yes	No	1~8	Every time you press [▲] or [▼], the value is increased or decreased by 1. The default value is: 6
Inhalation time	No	No	No	No	Yes	Yes	0.5~3.0 S	Every time you press [▲] or [▼], the value is increased or decreased by 0.1 S. The factory default is: 3.0S
Respiratory rate	No	No	No	No	Yes	Yes	6~40BPM	Every time you press [▲] or [▼], the value increases or decreases 1BPM. The default is: 10BPM

Backup ventilation	No	No	No	Yes	No	No	OFF/ ON	When the equipment is in S mode, the backup ventilation function is turned on and the patient has no spontaneous breathing within 6S, the backup ventilation frequency is 10 times per minute. The default is OFF
Automatic Opening	Yes	Yes	Yes	Yes	Yes	Yes	OFF/ ON	This function enables the device to automatically start and output preset pressure after the patient wears a breathing mask and takes several deep breaths. Press [▲] or [▼] to change settings. The default is: ON
Automatic shutdown	Yes	Yes	Yes	Yes	Yes	Yes	OFF/ ON	This function enables the device to automatically stop output pressure and shut down after the patient takes off the respiratory mask. Press [▲] or [▼] to change the settings. When leaving the factory, the default is: OFF
Tips for air leakage	Yes	Yes	Yes	Yes	Yes	Yes	OFF/ ON	When the patient's tube or mask falls off or breaks down during use, the device will give an alarm, prompting the patient to wear it again (this function is effective only when automatic shutdown is prohibited). Press [▲] or [▼] keys to change settings. The default is OFF
Expiratory decompression	No	Yes	No	No	No	No	OFF/ ON	During the use, the breathing rhythm is automatically detected, and the pressure in the mask is lowered when exhaling to make the patient feel comfortable. Press [▲] or [▼] to change the settings. When leaving the factory, the default is OFF

Notes: The following are the relationships among the parameters:

Minimum inspiratory pressure < inspiratory time < maximum inspiratory time

Initial pressure < expiratory pressure < minimum inspiratory pressure < inspiratory pressure < maximum inspiratory pressure.

3.Setting Options and Description of Service Menu

This menu is open only for authorized customer service engineers of the company (Code: 1709).

Setting an option	Set range	Description
Unit	hPa/cmH ₂ O	Used to select pressure measurement units, hPa or cm water column (cmH ₂ O). Press [▲] or [▼] to switch. The factory default is: hPa
Time format	12-hour/24-hour system	Two formats for selection, press [▲] or [▼] to switch, factory default is: 12 hours system
Date format	Year/month/day, Month/day/year, Day/month/year	There are three formats to choose from. Press [▲] or [▼] to switch. Factory default is: month / day / year
LCD backlight	Automatic/constant brightness	There are two modes of backlight for LCD display. Press [▲] or [▼] to switch. Factory default is: automatic
Running time	-----	Actual running cumulative time of display equipment
Pressure calibration	-----	Start-up pressure correction
Mask test	Start mask test	Test the air leakage of the mask.
Language	English/Chinese	There are two languages to choose from. Press [▲] or [▼] button to switch. Factory default is: English
Restore factory settings	-----	Press [▲] or [▼] key to choose whether to restore the default settings.
Clear User Data	-----	Press [▲] or [▼] key to choose whether to clear user data in SD card.

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Modèles: BR-20C, BR-20A, BR-20S, BR-20ST, BR-25C, BR-25A, BR-25S, BR-25ST

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Novidion GmbH
Shenzhen Bi-rich Medical Devices CO., LTD
SUNGO Cert GmbH
SUNGO Technical Service GmbH
Novidion GmbH, Fuggerstr. 30, 51149 Köln, Germany