DE Problembehebung

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EN Troubleshooting

ACP-25 Sleep apnea respiratory therapy device (Model: BR-25ST)

Problems related to the device:

Problem	Possible causes	Solutions
The mask pressure differs from the actual pressure	Improper hose connection	Connect the hose to the correct connection point.
	A leak in a mask or pressure sensor	Please contact the supplier of the equipment.
Devices do not work	Automatic switch-off function activated	Press the power button to restart the device.
	Incorrect insertion of the power plug	Insert the power plug correctly.
	Poor contact of the power plug	Please consult a professional electrician to have the outlet replaced.
Standby status cannot be started	Procedures for treatment devices were not started properly	Disconnect the device from the power supply, wait 30 minutes before switching the device back on, and then restart it.
Excessive noise	Improper hose connection	Connect the hose to the correct connection point.
The output temperature of the sensory device is abnormal	The air path of the unit may not be level and the incoming air may not be sufficient.	Replace air filter cotton and clean air vent
		Stay away from objects that may impede air circulation.
Low outlet air pressure	The air hose could be blocked	Replace the air filter cotton and clean the vent to check that the unit's intake path is clear.
	When the delayed booster function is used, it takes time for the pressure to rise to a certain level. Normal if the time is within the range.	If required, the delayed booster function can be switched off or a shorter delay of the boost time can be set.

Problems related to the patient:

Problem	Possible causes	Solutions
Dry mouth and throat	Non-heating humidifier	Switch on the heater and the humidifier.
	The temperature of the heater and humidifier is too low	Increase the temperature of the heating and humidification device
Red face	The mask fits too tight	The mask should fit as loosely as possible without allowing air to escape.
	The size of the mask is not correct	Every person's face is different. You can consider changing the mask.
	Allergy to mask materials	Consult a doctor or the supplier Place a skin-friendly material between the skin and the mask.
Shortness of breath	Excessive setting of treatment or expiratory pressure.	Consult a physician to appropriately reduce treatment pressure.

Difficulties with inhalation	Air leakage and hose rupture.	Contact the after-sales service personnel to replace the tube.
	Treatment or inhalation pressure set too low.	Consult the doctor and set the treatment pressure correctly.
	Inhalation sensitivity set too low.	Please choose a highly sensitive inhalation.
	The inlet is blocked.	Remove the object blocking the inlet.
Pain in the nose or ear	Inflammation of the paranasal sinuses or middle ear	Contact your doctor immediately
Unable to adapt to treatment pressures, feeling of discomfort	If the treatment pressure is higher than 13 hPa, patients feel uncomfortable. But the pressure is adjusted according to the condition and too low pressure cannot prevent apnea.	Patients need to get used to the compressed air for up to 4 weeks. They can try to relax and learn to breathe through their nose. If the problem persists, please contact your doctor.

Meaning of the warning buzzer:

Alert level	Sound description
High priority	dididi didi dididi didi
Medium priority	dididi dididi
Low priority	di di

Type	Warning priority	Description
Power Error!	high	Power Error: If the power supply is not normal, the microcontroller cannot operate properly and may not be able to send an alarm.
High temperature of motor!	high	Motor temperature: If the switching MOSFET is punctured, the temperature of the motor may be too high.
Leak!	medium	Leakage: A leakage exceeds 65 l/min and lasts longer than 36 seconds.
Low pressure!	medium	Low pressure: The treatment pressure is less than 3 hPa and lasts longer than 30 seconds.
SD card error!	medium	MicroSD card error: If the SD card cannot be written to, an alarm is triggered.

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