INHALATOR PULOX IN-200 VAPO

Instruction manual – Mesh Nebulizer pulox IN-200 Vapo

www.pulox.de

Inhalt

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1.3 Components and Accessories	30
2. Specifications Basic Specification.....	32
3. Safety Precaution.....	32
3.1 Warnings	33
3.2 Caution	33
4. Using Instructions.....	33
4.1 Connect the power supply	34
4.2 Fill the medication.....	36
4.3 Connect Medicine Cup.....	37
4.4 Install the Inhalation Mask.....	37
4.5 Wear the Mask.....	37
4.6 Start Inhalation	38
4.7 Stop inhalation	39
5. Cleaning and Disinfection.....	40
5.1 Clean the Medicine Cup and Mask	40

5.2 Disinfect the Medicine Cup and Mask	40
5.3 Replace the Medicine Cup and Mask.....	41
6. Maintenance	41
7. Trouble Shooting.....	42
8. Disposal	43
9. Electromagnetic Compatibility.....	43

Novidion GmbH
Fuggerstr. 30 • 51149 Köln
Tel.: +49 (0) 2203 / 9885 200
Fax: +49 (0) 2203 / 9885 206
www.pulox.de • Mail: info@novidion.de

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Operating Manual

Dear customer,
we are pleased that you have decided to purchase a product from our range.

Please read these operating instructions carefully and keep them accessible for other users. Pay attention to the safety precautions and notes. We are not responsible for any damage caused by failure to follow the instructions for use. We reserve the right to update and modify the instructions for use without prior notice.

With kind regards

Your Novidion-Team

1. Device Description

1.1 Product Instruction

Working principle: thanks to self-excited oscillation circuit, the ceramic sheets contained in nebulizing head together with the Metal Mesh are driven successively to oscillate at a high frequency, which makes the liquid flow through the micromesh of the Metal Mesh and be ejected and nebulized to form inhalant aerosol. Through Inhalation Mask to the patient's respiratory system, in order to achieve the purpose of inhalation treatment.


1.2 Symbols

The warning signs and the sample icons shown here are listed for you to use this product safely and correctly as well as to prevent the risk and the damage to you and others.

	Manufacturer		Operating instructions
	The device should be sent to the special agencies according to local regulations for separate collection after its useful life.		Type BF Applied Part
	Class II		Manufacturer's batch code

1.3 Components and Accessories

Before installing and debugging the product, the user first checks whether the appearance of the product is good, whether the variety and quantity of the accessories are consistent with the list of accessories after the instructions, and if there is any defect, please contact the supplier or manufacturer in time.

1	Medicine Cup	2	Mask
3	String	4	Battery Box (optional)
5	Main Unit	6	Main Unit Power Button
7	Main Unit Connection Cable	8	Battery (optional)
9	Type-C Adapter Head (optional)	10	Lightning Adapter Head (not included)
11	Micro-USB-Adapter Head (optional)		

2. Specifications

Basic Specification

Product Name	VAPO Mesh-Vernebler
Model	VP-M10/IN-200 Vapo
Power Supply	2X1,5V AA-alkaline batteries or DC 5V 1A
Nebulization Rate	≥0,20ml/min
Particle Size	MMAD approximately 5um
Sound Level	≤50dB
Ultrasonic Frequency	Approximately 114kHz
Capacity of Medicine Cup	≤8ml
Water Temperature of Medicine Cup	≤60°C
Weight of Medicine Cup	Approximately 11g
Dimensions	42(W)mm x 34(H)mm x 40 (D)mm
Safety Classification	Equipment that cannot be used in the case of flammable narcotic gas mixed with air or flammable narcotic gas mixed with oxygen or nitrous oxide
Electrical safety requirements	internal power supply equipment, BF application part
Operation mode	continuous operation
Working Environment	Temperature: 10 to 40°C Humidity: 30 to 85%RH Atmospheric Pressure: 86kPa~106kPa
Transport and Storage	Temperature: -20 to 50°C Humidity: 30 to 85%RH

Environment	Atmospheric Pressure: 86kPa~106kPa
-------------	------------------------------------

3. Safety Precaution

3.1 Warnings

- For this type, dose and regimen of the medication, be sure to follow the instructions of a doctor. Your symptom(s) may worsen.
- Clean and disinfect Medicine Cup and Inhalation Mask, after each use.
- Be sure to dry the cleaned and disinfected parts promptly, and store them in a clean place. Miscellaneous bacteria may propagate and you may get infected.
- This Unit is designed for single personal use. It is not appropriate for multi-users.
- Keep the Unit out of reach of children or infants. They may swallow small pieces. The Main Unit may be damaged or may cause injury to children or infants. In case a child swallows small pieces, immediately consult with a doctor.
- Do not wash the Main Unit and Battery Box with water or hot water. Do not immerse them in water. The Unit may experience electricity leakage or you may suffer electric shock. The Unit will cause trouble and may not be able to be used.
- This device is a professional atomization treatment equipment. Children and adults who are unable to take care of themselves should ensure that a guardian is present when using it.

3.2 Caution

- Do not inhale by using water in the Medicine Cup. Your symptom(s) may worsen.
- Do not drop medication on the Main Unit. If you drop medication, immediately wipe it off with gauze. If you use the Unit while it is wet, it may cause trouble or you may suffer electric shock.


- Do not poke Metal Mesh with a cotton swab or a pin. The Unit may be broken and cannot be used.
- Do not drop or apply strong shock to the Battery Box or Medicine Cup. They may be broken and cannot be used or you may suffer electric shock.
- Do not use the Battery Box other than the one supplied by VAPO. Do not use a broken Main Unit Connection Cable.
- In order to ensure the correct use of this product, please read this manual carefully and operate in strict accordance with the requirements of the manual. If you have any questions, please contact the supplier or manufacturer.
- When the storage and transportation temperature is lower than 5°C, the device should be restored to normal working conditions before use.

4. Using Instructions

4.1 Connect the power supply

The equipment supports external power supply (5V 1A) and internal battery box power supply (2X1.5V AA alkaline batteries). Only one of them can be selected for power supply. But the Battery Box is optional.


4.1.1 The Unit uses external power supply, such as Power Bank, Adapter. Be noted that this Main Unit is USB Connector and Power Supply is 5V 1A. Plug the USB connector of the Main Unit into the Power bank.


Power bank or Adapter is not included in the complete set of the product, which meets electrical grade and passes IEC 60601-1 certification.

4.1.2 The Unit uses external power supply with Battery Box. Be noted that rechargeable batteries cannot be used in the Battery Box!


1. Open the battery cover.

2. Insert batteries correctly by aligning the polarities (+ and -) with the battery indication mark.

3. Plug the Main Unit into the Battery Box and confirm the full connection tightly.


4.1.3 Battery life and battery replacement

- When using alkaline batteries: The Unit can be used for approximately 3 days (If it is used for 30 minutes a day).
- When the red light is on, batteries will be worn out soon. Replace both batteries with new ones.
- If you are not going to use the unit for a long period of time (3 months or more), remove the batteries.
- If you leave the Unit without using while keeping the batteries in it, the batteries may leak and cause trouble to the Unit.
- Do not insert batteries with the polarities (+ and -) in the wrong direction.

4.2 Fill the medication

Open the Medicine Cup, fill in the medication, close the Medicine Cup tightly. When using high concentration, high viscosity, suspension and volatile liquid, be sure to follow the doctor's instructions.

4.3 Connect Medicine Cup

Insert the other end of the Main Unit Connection Cable into the input port of the Medicine Cup and make sure the tight connection. Do not shake violently when the Medicine Cup is filled with liquid.

4.4 Install the Inhalation Mask


4.5 Wear the Mask

4.6 Start Inhalation

Short press the power button of the Main Unit to start inhalation and the light of the Main Unit is on.


- By taking a slow, deep breath inhale drugs.
- For the agents with high surface activity or high viscosity, the nebulization rate may be reduced.
- It is normal that a certain amount of sound will be produced by the Metal Mesh during nebulization.
- If the solution has high surface activity, a small amount of the solution may exudate the Metal Mesh.
- If the liquid in the Cup runs out, press the switch button on the Main Unit to turn off the device to avoid potential damage to the Metal Mesh.
- When the Medication is about to run out, users are advised to tilt the cup slightly to themselves to facilitate the remaining liquid to touch the Metal Mesh for nebulization.
- When inhaling with a Mask, do not cover the Mask vent with your hands or other objects.
- During aerosol inhalation. There is a vent on the lid of the Medicine Cup, please do not cover it with your hands or other objects, so as

not to affect the normal spray.

- During the inhalation, if the user feels uncomfortable, please turn off the machine immediately and consult the doctor in time.

4.7 Stop inhalation


Short press the power button of the Main Unit to stop inhalation.

5. Cleaning and Disinfection

5.1 Clean the Medicine Cup and Mask

- 1) Remove Inhalation Mask and Medicine Cup from the Main Unit.
- 2) Open the Medicine Cup and discard the medication.
- 3) Add a small amount of (2~5mL) water to the Medicine Cup, then fasten the lid tightly. Gently shake the Cup so that the residual liquid is fully dissolved in water. Open the Medicine Cup, pour the solution out of the Cup.
- 4) Wash the Mask in the water.
- 5) After washing the parts, please dry thoroughly in a clean environment.
- 6) Store all parts of the nebulizer in a clean place.

5.2 Disinfect the Medicine Cup and Mask

- 1) Use 75% medical alcohol to scrub the Mask, Medicine Cup for more than 3 minutes, then dry the Mask, Medicine Cup and storage them in the clean area..
- 2) After alcohol disinfection, rinse thoroughly with water to avoid alcohol residue. Be sure to disinfect the Medicine Cup and Mask after each used, and replace the parts in time if the parts are seriously contaminated.

5.3 Replace the Medicine Cup and Mask

The Medicine Cup should be replaced every 3 months at least.

Otherwise, functional capacity and safe treatment will no longer be guaranteed. It also need to be replaced with a new Medicine Cup if the Medicine Cup is damaged or cannot work normally.

6. Maintenance

- Before each time used, please check the device carefully. Do not use the device if it is damaged in any way. The continuous use of a damaged unit may cause injury, improper results, or serious danger.
- Store and use the device at a cool, dry and ventilated environment. Avoid to approach to the fire and the heat sources, or it will cause the battery to explode.
- If you have any problems with this device, such as setting up, maintaining or using, please contact with SERVICE PERSONNEL of Vapo Healthcare Co., Ltd. Don't open or repair the device by yourself.
- Please report to Vapo Healthcare Co., Ltd. if any unexpected operation or events occur.

7. Trouble Shooting

Trouble	Possible Cause	How to correct
The nebulization rate is extremely low.	The Metal Mesh is stained.	Clean the Medicine Cup to remove the stain. If the trouble cannot be solved by the above-mentioned procedure, replace the Medicine Cup with a new one.
	The red light of the Main Unit is on. This shows that the battery voltage is low and the battery is exhausted.	Replace the batteries with new ones.
The Unit does not nebulize.	The Main Unit is not correctly connected.	Remove the Main Unit from the Battery Box and connect it correctly.
	The Main Unit power is not turn on.	Turn on the Main Unit power switch.
	The polarities (+ and -) of batteries are in the wrong direction.	Insert batteries in the correct alignment.
	The batteries with low power.	Replace batteries with new ones.
The light of main unit works correctly but the Unit does not nebulize.	Main Unit connection cable damaged.	Replace the new Main Unit.
	The Metal Mesh is extremely stained.	Clean the Medicine Cup according to the

		instructions. If there is no troubleshooting after cleaned, please replace the Medicine Cup with a new one.
Liquid leakage.	The Medicine Cup is damaged.	Replace a new Medicine Cup.

8. Disposal

DO NOT dispose the batteries in domestic waste.

Dispose of the batteries according to the local regulations dealing with the disposal of these special materials (e. g. to the collecting points).

The device is made of a combination of plastic and stainless steel. Be sure to dispose of it in accordance with local regulations as unsorted municipal waste. You may recycle it at your local community or appliance recycling center.

9. Electromagnetic Compatibility

Guidance and Manufacturer's Declaration

1) This equipment needs to be installed and put into service in accordance with the information provided in the ACCOMPANYING DOCUMENTS;

2) IMMUNITY TEST LEVELS for basic safety and essential performance of ME equipment and ME systems should be chosen based on a high

probability of maintaining basic safety and essential performance, and shall be according to the professional Healthcare facility environment, home Healthcare environment, and special environment, based on the locations of intended use.

3) HOME HEALTHCARE ENVIRONMENT is dwelling place in which a patient lives or other places where patients are present, excluding professional Healthcare facility environments where operators with medical training are continually available when patients are present. Such as schools, outdoors, domiciles, vehicles hotels and pensions.

EXAMPLE: As indicated in Table 6 of IEC 60601-1-2: 2014 for ME EQUIPMENT, a typical cell phone with a maximum output power of 2W yields $d = 3.3\text{m}$ at an IMMUNITY LEVEL of 3V/m..

A1 Electromagnetic Emissions - For all EQUIPMENT and SYSTEMS

Guidance and manufacture's declaration-electromagnetic emission		
The VP-M10/IN-200 Vapo is intended for use in the electromagnetic environment specified below; The customer or the user of the VP-M10/IN-200 Vapo should assure that it is used in such and environment.		
Emission test	Compliance	Electromagnetic environment-guidance
RF emissions CISPR 11	Group 1	The VP-M10/IN-200 Vapo uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to

		cause any interference in nearby electronic equipment.
RF emissions CISPR 11	ClassB	The VP-M10/IN-200 Vapo is suitable for use in all establishments, other than domestic and those directly connected to the public low-voltage power supply network that supplies building used for domestic purposes.
Harmonic emissions IEC 61000-3-2	ClassA	
Voltage fluctuations/-flicker emissions IEC 61000-3-3	Complies	

A2 Electromagnetic Immunity - For Home Healthcare Environment EQUIPMENT and SYSTEMS

Guidance and manufacture's declaration - electromagnetic immunity			
The VP-M10/IN-200 Vapo is intended for use in the electromagnetic environment specified below; The customer or the user of the VP-M10/IN-200 Vapo should assure that it is used in such and environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment -guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8kV contact ±2kV, ±4 kV, ±8kV, ±15kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	Floors should be wood, concrete or ceramic tile. If floor are covered with synthetic material, the relative humidity

			should be at least 30%.
Electrical fast transient/ burst IEC 61000-4-4	±2kV for power supply lines	±2kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1kV Leitung zu Leitung ±2kV Leitung zu Erde	±1kV Leitung zu Leitung ±2kV Leitung gegen Erde	Mains power quality should be that of a typical commercial or hospital environment.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	30A/m	30A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % UT; 0,5 Zyklus Bei 0°, 45°, 90°, 135°, 180°, 225°, 270° und 315°	0 % UT; 0,5 Zyklus Bei 0°, 45°, 90°, 135°, 180°, 225°, 270° und 315°	Mains power quality should be that of a typical commercial or hospital environment. If the user of the VP-M10 requires continued operation during power mains interruptions, it is recommended that the VP-M10 be powered from an uninterruptible power supply or a battery.

<p>ConductedHF IEC 61000-4-6</p> <p>RadiatedHF IEC 61000-4-3</p>	<p>3Vrms 150kHz bis 80MHz</p> <p>10V/m 80MHz bis 2,7GHz</p>	<p>3Vrms</p> <p>10V/m</p>	<p>The MANUFACTURER should consider reducing the minimum separation distance, based on RISK MANAGEMENT, and using higher IMMUNITY TEST LEVELS that are appropriate for the reduced minimum separation distance. Minimum separation distances for higher IMMUNITY TEST LEVELS shall be calculated using the following equation:</p> $E = \frac{6}{d} \sqrt{P}$ <p>Where P is the maximum power in W, d is the minimum separation distance in m, and E is the IMMUNITY TEST LEVEL in V/m.</p> <p>NOTE UT is the a.c. mains voltage prior to application of the test level.</p>
--------------------------------------------------------------------------	-------------------------------------------------------------------------------------	---------------------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Vertrieb durch/ Distribution by/ Distribución por/ Distribué par/
Distribuito da:

Novidion GmbH, Köln

Hersteller/Manufacturer/Fabricante/Producteur/Produttore:

Vapo Healthcare Co., Ltd.

Southern unit of third floor, building B, No. 99 Yudai West Rd,

High tech district, KunShan, Suzhou, Jiangsu Province, China

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Share Info GmbH
Novidion GmbH, Fuggerstr. 30, 51149 Köln