SONOTRAX C ULTRASCHALL
FETAL DOPPLER - BEDIE-
NUNGSANLEITUNG

Sonotrax Ultrasound Fetal Doppler
- Instruction manual

DE
Novidion GmbH
Fuggerstr. 30 • 51149 Köln • Deutschland
Tel.: +49 (0) 2203 / 9885 200
Fax: +49 (0) 2203 / 9885 206
www.pulox.de • Mail: info@pulox.de

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Instruction manual

Dear customer,

We are pleased that you have chosen a product from our range. Please read these instructions for use carefully and keep them accessible to other users. Please follow the instructions.

With kind recommendation

Your Novidion team

Note for users

Read these operating instructions carefully and thoroughly and follow the instructions for use exactly. Failure to follow the operating instructions may result in measurement errors, damage to the device and personal injury. The manufacturer is not responsible for safety, reliability and performance problems and any other monitoring deviation, personal injury or equipment damage due to user negligence. The manufacturer's warranty does not cover this.

This product is a medical device that can be used repeatedly.

Note: Read these operating instructions carefully before using the product.

In the event of changes and software updates, the information in this manual may be changed without prior notice.

Attention

This user manual has been written and prepared in accordance with the Medical Device Directive and harmonized standards (MDD93/42/EEC). In the event of changes and software upgrades, the information in this document is subject to change without notice.

The manufacturer makes no warranty of any kind with respect to this material, including, but not limited to, the implied warranties of merchantability and fitness for a particular purpose. The manufacturer assumes no responsibility for any errors in this document or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

No part of this document may be copied, reproduced or translated into another language in any form without the express written permission of the manufacturer.

The information in this document is subject to change without notice.

Responsibility of the manufacturer

The manufacturer only considers itself responsible for the safety, reliability and function of the devices if:

Installation and repairs are carried out by persons authorized

by the manufacturer and the appliance is used in accordance with the instructions for use.

 WARNING

This device is not intended for treatment. It is used to record the fetal heart rate (FHR). If there is any doubt about the fetal heart rate result, please use other methods for immediate verification, e.g. a stethoscope.

Warranty

The device must not be repaired by the user. All maintenance work must be carried out by qualified personnel authorized by the manufacturer. We warrant that every product sold is free from defects in workmanship and materials and conforms to the product specifications stated in this manual. If the product does not function as warranted within the warranty period, we will repair or replace the device free of charge. Improper use or improper maintenance will void the warranty and guarantee.

Use of notes

Notes indicate the most important precautionary measures.

 **WARNING:** The WARNING label warns of certain actions or situations that could lead to personal injury or death.

 **CAUTION:** The CAUTION label warns of actions or situations that may damage the device, result in inaccurate data or invalidate a procedure.

Note: A NOTE provides useful information about a function or procedure.

Chapter 1: Security guidelines

The device belongs to the internally operated device type; the applied part has a shock protection rating of type CF.

Type B application protection means that the patient connections comply with the permitted leakage currents and the electrical voltages comply with IEC 60601-1.

1.1 Safety precautions

WARNING and CAUTION instructions must be observed. To avoid the risk of injury, please observe the following precautions when using the appliance.

 Warning notices:

- This device is not explosion-proof and must not be used in the vicinity of flammable anesthetics.
- Do not throw batteries into a fire as they could explode.
- Do not attempt to charge normal dry batteries as they may leak and cause a fire or explode.
- Do not touch the input or output connection and the patient at the same time.
- The device accessories that are connected to the analog or digital interface must be certified in accordance with the relevant IEC standards (e.g. IEC 950 for data processing accessories or IEC 60601-1 for medical

accessories). Furthermore, all configurations must comply with the valid version of the IEC 60601-1-1 system standard. Anyone connecting additional accessories to the signal input or output connector is configuring a medical system and is therefore responsible for ensuring that the system complies with the requirements of the applicable version of the IEC 60601-1-1 system standard. If in doubt, please contact our technical service department or your local distributor.

- Fetal Doppler ultrasound is a tool for healthcare professionals and should not be used in place of normal fetal monitoring.
- The batteries should only be changed outside the patient's environment (1.5 meters away from the patient).
- Please use the ultrasound fetal Doppler probe provided by the manufacturer.
- Do not pull the cable of the probe further than 2 meters, otherwise the probe could break out of the connection socket of the ultrasound fetal doppler.
- The appliance may only be serviced by authorized and qualified personnel.

CAUTION:

- The device has been designed for standard continuous operation. Do not immerse it in liquids (not drip or splash-proof).

- Protect the device from dirt. Avoid shocks.
- Do not use high-temperature methods, electron beam or gamma ray disinfection for disinfection.
- Electromagnetic interference - Make sure that there are no sources of strong electromagnetic interference, e.g. radio transmitters, cell phones, etc., in the area where you are using the device.
- Before use, the user must check that the device has no visible faults that could affect patient safety or monitoring performance. The recommended inspection interval is one month or less. If there is an obvious defect, replacement of the device is recommended.
- The following safety checks should be performed once every 2 years or as specified in a test and institution inspection protocol by a qualified person who has adequate training, knowledge and practical experience to perform tests.
 - * Check the device for mechanical and functional damage.
 - * Check the safety-relevant markings for legibility.
 - * Check that the appliance is working properly as described in the operating instructions.
 - * Testing the patient leakage current in accordance with IEC 60601-1/1988: Limit: μA (B).

The leakage current should never exceed the limit. The data should be recorded in an operating log. If the appliance does not function properly or fails any of the above tests, the appliance must be repaired.

CAUTION

- The battery must be disposed of properly after use in accordance with local regulations.
- The battery must be removed from the battery compartment if the device is not to be used for a longer period of time.
- Store the batteries in a cool and dry place.
- To ensure the performance and service life of rechargeable batteries, please fully charge the batteries before using them for the first time. Normally, batteries must be charged continuously for 14 hours or according to the instructions on the batteries.
- Do not insert the anode and cathode of the batteries the wrong way round.
- The specified period of use of this product is 5 years.
- Please return the device to the manufacturer at the end of its service life or dispose of it in accordance with local regulations.
- The device must not be used together with a defibrillator or high-frequency surgical equipment.

- Please only use accessories that have been authorized by our company or the device could be damaged.
- Please keep the probe away from sharp-edged tools.
- Please use the ultrasound Fetal Doppler in an environment that is free of strong electromagnetic fields that could influence the measurement results.

Cleaning the appliance:

ACHTUNG:

- Do not use strong solvents, e.g. acetone.
- Never use abrasives such as steel wool or metal polish.
- Do not allow any liquids to enter the appliance. Do not immerse any parts of the appliance in liquids.
- Avoid pouring liquids onto the appliance during cleaning.
- Do not leave any cleaning agents on the surface of the appliance.

Disinfection of the appliance:

WARNING:

- Never disinfect the probe or accessories using low-temperature steam disinfection or other methods (observe the safety instructions in the accompanying documents).

Chapter 2 Introduction

2.1 Overview

The Ultrasound Fetal Doppler is a portable obstetric device that can be used in the hospital, clinic or at home for daily self-monitoring by the pregnant woman.

The device has a colorful high-resolution LCD screen that displays the fetal heartbeat waveform and heart rate to assist the doctor's timely diagnosis.

It includes components of ultrasonic signal transducer and receiver, analog signal processing, FHR calculation, LCD screen control, etc.

It has 3 operation modes: real-time FHR display mode, average FHR display mode and manual mode.

It has an audio output and can be connected to headphones or recording devices with audio input.

2.2 Properties

- * Portable Fetal Doppler
- * Curved probe for easy use
- * Colorful screen with fetal heart rate value, bar graph and heart rate curve
- * Fetal heart rate is displayed in red if the value is not within the normal range
- * Battery level indicator
- * Probe can be changed
- * Probe check

- * Built-in loudspeaker
- * Output for headphones
- * Automatic switch-off function
- * 2 standard 1.5 V alkaline batteries

Chapter 3 Overview and configuration


Fig. 3-2 Rear panel



Fig. 3-3 Top side

1. Loudspeaker	1. Headphones
2. Battery compartment	2. Headphone socket
	3. Caution

1. Display
2. Mode Selection
3. Power Button
4. Probe
5. Start/Stop

Fig. 3-1 Front side

3.1 Display

The LCD display of the ultrasound fetal doppler:

1. FHR value
2. Mode
3. Nominal frequency of the probe
4. Battery level indicator
5. FHR bar display
6. FHR heartbeat curve


Fig. 3-4 LCD display

3.2 Operating buttons

There are 3 buttons (power, mode and start/stop) and a volume control. These have the following functions:

3.2.1 Power button

Function: Switching on and off

On: Press the button once

Off: Press the button and hold it down for approx. 3 seconds.

3.2.2 Mode button

Function: Selecting the mode

Press once to switch to the next mode during use.

The ultrasound fetal doppler has a reminder function. When you switch the device on, it automatically starts in the last mode used.

3.2.3 Start/Stop/Taste

Function: If you press the button in mode 3, the fetal heart rate starts to count. Press the button

again to stop counting.

3.2.4 Volume control

Rotary control for adjusting the volume.

Turn the control from left to right to adjust the volume from loud to quiet.

3.3 Description of the top of the device

Headphone connection: The audio output can be connected to headphones or a recording device with an audio input.

: The socket or connection to which the plug of the headphones is connected.

: Attention. Please observe the enclosed documents.

The accessories that are connected to the analog or digital interface must be certified according to the applicable IEC standards (e.g. IEC 950 for data processing accessories and IEC 60101-1 for medical accessories). Furthermore, all configurations should comply with the valid version of the IEC 60601-1-1 system standard. Anyone who connects additional accessories to the signal input or output configures a medical system and is therefore responsible for ensuring that the system complies with the conditions of the valid version of the IEC 60101-1-1 system standard. If in doubt, please contact our technical service department or your local distributor.

Signal interface


Fig. 3-5 Headphone connection for audio output

Headphone connection, as shown in Fig. 3-5, the description of the plugs is as follows:

Stecker	Beschreibung
1	GND

2	Signal
3	Signal
4	Signal
5	Signal

3.4 Ultrasound probe

2MHz Ultrasound probe: The 2 MHz ultrasound probe is the standard probe supplied for fetal Doppler.


Fig.3-6 2MHz Ultrasonic probe

The CD2.0 marking on the probe has the following meaning:

C: Operating mode of the probe is Continuous Wave (CW)

D: The structural form of the probe is cell type.

2.0: The frequency of the probe is 2 MHz.

Chapter 4 General operation

4.1 Fetal heart rate monitoring

1. Switch on the device by pressing the power button (the LCD display is as shown in Fig. 3-4).
2. Find the position of the fetus: first feel the position of the fetus with your hand. Find the best orientation to monitor the fetus. Apply a generous amount of gel to the probe and place the probe where you can best monitor the fetal heart. Adjust the probe by moving it back and forth so that you receive an optimal audio signal. Adjust the volume according to your requirements.
3. Fetal heart rate calculation: The LCD screen will now display the fetal heart rate as a value, the bar graph and the heartbeat curve.
4. Switch off the device: Press and hold the power button for approx. 3 seconds to switch off the device.

CAUTION:

1. Place the probe in the position where you can best recognize the heart rate. This will give you better results.
2. Do not place the probe where there is a strong placental blood flow sound (PBS) or a strong umbilical sound (UMS).
3. If the pregnant woman is in a horizontal position and the fetus is in a normal position, place the probe at the

midline below the umbilicus to obtain the clearest FHR sound.

4. Do not record the fetal heart rate unless an audible fetal sound is heard.

4.2 Mode selection

4.2.1 Real-time display mode of the fetal heart rate (mode 1)

While you are recording the fetal heart rate signal, the bar graph on the LCD screen shows the strength of the fetal heartbeat signal. At the same time, the real-time heart rate is displayed as a value and the fetal heartbeat as a curve.

4.2.2 Average display of the fetal heart rate (mode 2)

This mode is used to obtain more stable heart rate readings. In this mode, the fetal heart rate averages 8 beats. As you capture the fetal heart rate signal, the bar graph on the LCD screen shows the strength of the fetal heartbeat signal. The displayed heart rate as a value and the fetal heart rate as a curve change slowly.

4.2.3 Manual mode

When you press the Start/Stop button in mode 3, the system starts counting the audible beats and the fetal heart rate is displayed in "-- -- --" format. Once the fetal heartbeat has been detected, the bar graph will display the strength of the heartbeat signal. Press the Start/Stop button again to finish the calculation and the device will automatically calculate the

fetal heart rate derived from the calculation time and display the result. This frequency value is retained until a new measurement is taken or the mode is changed.

4.3 Application of the probe

4.3.1 Checking the probe

If the probe slips out of the Fetal Doppler, the LCD display shows a flashing "-- -- -- --" and "Probe fall!". The frequency data display disappears. The probe must now be reconnected. As soon as the probe is reconnected, the "Probe fall!" displayed on the screen disappears and the frequency data of the probe is displayed again.

4.3.2 Replacing the probe

A probe has already been connected to the Fetal Doppler during packaging by the manufacturer. If you need to replace this with another probe, first switch off the ultrasound Fetal Doppler and then remove the probe from the holder of the Fetal Doppler. Then pull the probe plug out of the socket. Then connect the plug of the replacement probe to the socket.

Note: Place the currently unused probe down carefully and avoid dropping, knocks, etc. If the ultrasound Fetal Doppler is not used for a longer period of time, we recommend removing the probe plug from the socket of the Fetal Doppler and clamping the probe in the holder. Then pack the ultrasound fetal doppler with the probe into the packaging.

4.3.3 Removing and placing the probe

1. Remove the probe: Hold the device firmly with one hand and grasp the handle of the probe with the other hand to remove it.


Fig. 4-1 Remove the probe

2. Placing the probe: This works in reverse to removing the probe. Hold the device with one hand and the top of the probe with the other. Then press the probe into the holder.

4.4 FHR Alarm

The normal range of the fetal heart rate is between 120 BPM - 160 BPM, the LCD screen displays the heart rate as a numerical value in green. If the heart rate is too fast or too slow, or outside the normal range, the numerical value will

turn red. To be on the safe side, have yourself checked by a doctor to ensure the safety of the fetus.

4.5 Battery status indicator

The LCD screen displays the battery status as follows:

	Battery is full
	Battery is not completely full
	Battery is low, replace the battery.

If the battery charge is too low for normal operation, the LCD screen will display "Low Power!". The symbol starts to flash and the appliance switches off automatically after a while.

4.5 Replacing the battery

The back of the device is at the top. First open the battery compartment and then remove the batteries from the battery compartment (see Figure 4-2).


Fig. 4-2 Change batteries

Insert two AA batteries into the battery compartment (follow the instructions in the battery compartment to insert the batteries correctly), then close the battery compartment again.

CAUTION:

The battery should be removed if the device is not to be used for an extended period of time.

Chapter 5 Explanation of symbols

Symbol	Description
	Type CF
	Observe operating instructions
	Headphone socket
	Volume control
	The CE marking certifies conformity with the essential requirements of Directive 93/42/EEC for medical devices
	Disposal in accordance with EC Directive 2002/96/EC - WEEE (Waste Electrical and Electronic Equipment)
	European representation
	Serial number

Chapter 6 Product specification

Product name: Ultrasound Fetal Doppler

Model: Sonotrax C

Safety: According to IEC 60601-1:1988+A1:1991+A2:1995

Classification:

Anti-electric shock type: Internally operated device

Anti-electric shock grade: Type CF Application part

Harmful liquids Degree of protection: Conventional device
(sealed device without liquid protection)

Degree of protection of the appliance in the presence of
flammable gases: Device not suitable for use in the presence
of flammable gases

Operating system: Continuous operation

EMC: Group I Class B

Recommended area of application: Suitable for use after the
12th week of pregnancy

Physical properties:

Dimensions: 135mm (L) x 92mm (W) x 35mm (H)

Weight: approx. 245 g (incl. batteries)

Environment:

Operation:

Temperature: +5°C ~ +40°

Humidity: ≤ 80%

Atmospheric pressure: 70kPa ~ 106kPa

Transportation and storage:

Temperature: $-10^{\circ}\text{C} \sim +55^{\circ}\text{C}$

Humidity: $\leq 93\%$

Atmospheric pressure: $50\text{kPa} \sim 106\text{kPa}$

FHR performance (fetal heart rate)

FHR measuring range: $50 \sim 240 \text{ BPM}$ (BPM: beats per minute)

Frequency resolution: 1 BPM

Accuracy: $\pm 2 \text{ BPM}$

Energy consumption: $< 0,8 \text{ W}$

Automatic switch-off function: After 1 min. without a signal.

Recommended battery type: 2 1.5 V DC batteries (size AA LR6)

2 MHz ultrasonic probe (standard equipment):

Nominal frequency: 2.0 MHz

Operating frequency: $2.0 \text{ MHz} \pm 10\%$

P-: $< 0.5 \text{ MPa}$

I_{ob}: $< 10 \text{ mW/cm}^2$

I_{spta}: $< 50 \text{ mW/cm}^2$

Ultrasonic output intensity: I_{sata} $< 5 \text{ mW/cm}^2$

Operating mode: Continuous Wave Doppler (CW Doppler)

Effective radiation area of the transducer: $208 \text{ mm}^2 \pm 15\%$

Note: In all operating modes mechanical index: $MI < 1$, thermal index: $TI < 1$

Chapter 7 Care and maintenance

7.1 Maintenance

The acoustic surface of the probe is fragile and must be handled with care.

The gel must be removed from the probe after each use.

These precautions will prolong the life of the device.

Before use, the user must check that the device shows no visible evidence of damage that could affect the safety of the patient or the performance of the Fetal Doppler. The recommended inspection period is one month or less. If damage is found, we recommend replacement before use.

The device should be subjected to regular safety tests to ensure proper patient isolation from electrical leakage currents. This should include measuring the leakage currents. The recommended test period is once every 2 years or as specified in the test or inspection protocol.

The accuracy of the FHR is controlled by the device and cannot be adjusted by the user. If the FHR result is not accurate, please use other methods, such as a stethoscope, to check immediately or contact your local dealer or the manufacturer for assistance.

7.2 Cleaning

Switch off the device and remove the batteries before cleaning.

Keep the outer surface of the device clean and free of dust and dirt, clean the outer surface (including the display screen) of the housing with a dry, soft cloth. If necessary, clean the housing with a soft cloth moistened with soap or water and wipe it dry immediately with a clean cloth. Clean the probe with a soft cloth to remove all residues of the ultrasonic coupling gel. Clean with soap and water only.

**CAUTION:**

- Do not use strong solvents such as acetone.
- Do not use abrasive cleaners such as steel wool or metal polish.
- Do not allow any liquid to run into the product and do not immerse any parts of the appliance in liquids.
- Avoid pouring liquids onto the appliance during cleaning.
- Do not leave any residues of cleaning agents on the surface of the appliance.

Notes:

Wipe the surface of the probe with 70% ethanol and allow it to air dry or clean it with a clean, dry cloth.

7.3 Disinfection

Clean the device housing, the probe, etc. as described above and then wipe the probe with an alcohol-impregnated wipe (70% ethanol).

Wipe the probe with a clean, dry cloth to remove any residual moisture.

WARNING: Never disinfect the probe or the device using low vapor or other methods.

Chapter 8 Solution to possible problems

If the following problems occur when using the ultrasound fetal doppler, solve them as follows:

Problems	Possible causes	Solution
Weak sound	Volume is too low Low battery level Gel has not been applied	Turn up the volume Change the batteries Apply gel
Noise	The probe is too close to the device Interference from external signals Low battery level	Increase the distance between the probe and the device Keep the device far enough away from external signals Change the batteries
Low sensitivity	Position of the probe is not correct Gel was not applied	Align the position of the probe Apply gel

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Anhang 3

Guidance and manufacture's declaration – electromagnetic emissions- for all EQUIPMENT and SYSTEMS

Guidance and manufacture's declaration - electromagnetic emission

The *Sonotrax C Pocket Fetal Doppler* is intended for use in the electromagnetic environment specified below. The customer of the user of the *Sonotrax C Pocket Fetal Doppler* should assure that it is used in such an environment.

Emission test	Compliance	Electromagnetic environment – guidance
---------------	------------	--

RF emissions CISPR 11	Group 1	The <i>Sonotrax C Pocket Fetal Doppler</i> uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emission CISPR 11	Class B	The <i>Sonotrax C Pocket Fetal Doppler</i> is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable	

**Guidance and manufacture's declaration – electromagnetic immunity –
for all EQUIPMENT and SYSTEMS**

Guidance and manufacture's declaration – electromagnetic immunity			
The <i>Sonotrax C Pocket Fetal Doppler</i> is intended for use in the electromagnetic environment specified below. The customer or the user of <i>Sonotrax C Pocket Fetal Doppler</i> should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±6 kV air	Floors should be wood, concrete or ceramic tile. If floor are covered with synthetic material, the relative humidity should be at least 30%. the manufacturer may recommend the ESD precautionary procedures to user.

Power frequency (50Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
--	-------	-------	---

Anhang 4

Allgemeine Empfindlichkeit													
Durchmesser des Zielreflektors (mm)	Abstand (d) (mm)	Reflexions-verlust A (d)	Zweiwegdämpfung $B=\sum B_a + B_a$							V_s (r.m.s) mV	V_n (r.m.s) mV	$C=20\log 10 \frac{V_s(r.m.s)}{V_n(r.m.s)}$	Allgemeine Empfindlichkeit (S=A(d)+B+C) dB
			$\sum B_a$ (T: mm B_a : dB)					B_a (dB)	B (dB)				
1,58 A=45,7 dB@2 MHz	50	45,7	T	20	4,8	4,0	-	0	57,6	186	94	5,93	109,2
			B_a	40	9,6	8,0	-						
	75	45,7	T	20	4,8	3,4	-	0	56,4	175	90	5,78	107,8
			B_a	40	9,6	6,8	-						
	100	45,7	T	20	4,8	3,4	-	0	56,4	174	89	5,82	107,9
			B_a	40	9,6	6,8	-						
	200	45,7	T	20	4,8	-	-	0	49,6	173	90	5,68	100,9
			B_a	40	9,6	-	-						
2,38 A=43,2 dB@2 MHz	50	43,2	T	20	4,8	3,4	2,2	0	60,8	178	89	6,02	110,0
			B_a	40	9,6	6,8	4,4						
	75	43,2	T	20	4,8	3,4	1	0	58,4	170	90	5,52	107,1
			B_a	40	9,6	6,8	2						
	100	43,2	T	20	4,8	3,4	-	0	56,4	165	85	5,76	105,3
			B_a	40	9,6	6,8	-						
	200	43,2	T	20	4,8	1	-	0	51,6	160	85	5,49	100,2
			B_a	40	9,6	2	-						
Doppler Frequenz (Hz)		333							Geschwindigkeit des Ziels (cm/s)		12,5		