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Instruction manual

Dear customer,

we are pleased that you have chosen a product from our range. Please read these operating instructions carefully and keep them accessible for other users. Observe the notes.

With kind regards

Your Novidion Team

Foreword

Please read the User Manual carefully before using this product. The User Manual which describes the operating procedures should be followed strictly. This manual describes in detail the individual steps that must be followed when using the unit, otherwise there is a risk of damage to the unit or injury to persons. Any anomalies or personal injury and device damage arising from use, maintain, store do not follow requirements of the User Manual, our company is not responsible for the safety, reliability and performance guarantees! The manufacturer's warranty service does not cover such faults! Our company has a factory record and user profile for each device, users enjoy free maintenance services for one year from the date of purchase. To enable us to provide you with a comprehensive and efficient maintenance service, we ask that you return the warranty card when repairs are carried out.

Note: Please read the User Manual carefully before using this product.

Described in this User Manual is in accordance with practical situation of the product. In case of modifications and software upgrades, the information contained in this document is subject to change without notice.

The warning items

Before using this product, you should consider the safety and efficacy of the following described:

- The individual measurement results are described in combination with clinical symptoms by qualified doctors.
- The reliability and operation of this product as described in this manual depends on following the maintenance instructions.
- The intended operator of this product may be the patient.
- Do not perform maintenance and service while the device is in use.

 **Warning:** Replace accessories which not provided by our company may lead to the occurrence of errors. Arbitrary replacement of adapters and cuffs can lead to incorrect measurement results. You should not attempt to service the product without trained service personnel from our company or other approved maintenance organisations.

Responsibility of operator

- The operator must read the user manual carefully before using the unit and strictly follow the operating instructions in the user manual.
- Security requirements must be fully considered during product use, but the operator should not neglect monitoring the patient and the condition of the machine.
- It is the responsibility of the operator to inform our company of the condition of use of the product.

Responsibility for our company

- Our company is responsible for delivering a qualified product that meets the company standard for that product.
- Our company will provide the circuit diagram, calibration method and other information at the request of the user to help the appropriate and qualified technicians to repair those parts designated by our company.
- Our company has the responsibility to complete product maintenance according to the contract.
- Our company has the responsibility to respond the requirements of user in time.
- In the following case, our company is responsible for the impact on the safety, reliability and performance of the device:
 - Assembly, addition, debugging, modification or repair are carried out by personnel approved by our company.
 - The electrical facilities in the room are in compliance with the relevant requirements and the device is used in accordance with the User Manual.

The User Manual is written by our company. All rights reserved.

Chapter 1 Functions and Purpose

1.1 Main Functions

- Measure blood pressure and store the measurement results.
- Data storage function, up to 199 records can be stored.
- With data review interface which is convenient for reviewing blood pressure parameter.
- The screen will prompt message when the power is low.

- When the measurement result cannot be obtained due to some factors during the measurement, the device will display the corresponding error information.
- Measurement units: mmHg and kPa, which can be switched by the button.
- With automatic shutdown function: If there is no operation, the device will automatically turn off.

1.2 Purpose

The Sphygmomanometer is intended to measure the non-invasive blood pressure of adult. It is designed to use in clinics, medical centers and daily blood pressure monitoring at home.

Chapter 2 Safety Precautions

In order to use it correctly, please read the "Safety Precautions" carefully before using it.

Operators do not need professional training, but should use this product after fully understanding the requirements in this manual.

To prevent users from suffering damage or loss due to improper use, please refer to "Safety Precautions" and use this product properly.

For safety reasons, be sure to comply with safety precautions.

Note

In case of improper use, there is a possibility of damage to persons and goods.

Good damage means the damage of house, property, domestic animal and pet.

Contraindication

No contradictions.

Warning

- You must not perform NIBP measurements in patients with sickle cell disease or in conditions where the skin is damaged or is likely to be damaged.
- In patients with severe coagulation disorders, automatic blood pressure measurement should be based on clinical assessment, as limb friction against the cuff may cause a risk of hematoma.
- For severe blood circulation disorder or arrhythmia patients, please use the device under the guidance of a doctor. If the arm is squeezed during measurement, it may cause acute internal hemorrhage or inaccurate measurement results.

Measurement Limitations

To different patient conditions, the oscillometric measurement has certain limitations. The measurement is in search of regular arterial pressure pulse. In

those circumstances when the patient's condition makes it difficult to detect, the measurement becomes unreliable and measuring time increases. The user should be aware that the following conditions could interfere with the measurement, making the measurement unreliable or longer to derive. In some cases, the patient's condition will make a measurement impossible.

Patient Movement

Measurements will be unreliable or cannot be performed if the patient is moving, shivering or having convulsions. These motions may interfere with the detection of the arterial pressure pulses. In addition, the measurement time will be prolonged.

Cardiac Arrhythmia's

Measurements will be unreliable and may not be possible if the patient's cardiac arrhythmia has caused an irregular heartbeat. This prolongs the measuring time.

Heart-lung Machine

Measurements will not be possible if the patient is connected to a heart-lung machine.

Pressure Changes

Measurements will be unreliable and may not be possible if the patient's blood pressure is changing rapidly over the period of time during which the arterial pressure pulses are being analysed to obtain the measurement.

Severe Shock

If the patient is in severe shock or hypothermia, measurements will be unreliable since reduced blood flow to the peripheries will cause reduced pulsation of the arteries.

Heart Rate Extremes

Measurements cannot be made at a heart rate of less than 40 bpm and greater than 240 bpm.

Obese patient

The thick fat layer of the body reduces the measurement accuracy, as the impulses cannot reach the cuff due to the damping.

Warning

- **Self-diagnosis and treatment based on measurement results can be dangerous.** Follow the instructions of your physician. Please hand measurement results to the doctor who knows your health and accept diagnosis.
- **For Infant and the person who can't express oneself, please use the device under the guidance of a doctor.** Otherwise, it may cause accident or dissension.
- **Please do not use for any other purpose except BP measurement.**

Otherwise, it may cause accident or holdback

- **Please use special cuff.** Otherwise, it is possible that the measurement result is incorrect.
- **Please do not keep the cuff in the over-inflated state for a long time.** Otherwise, it may cause risk.
- **Do not use the device if there are flammable anesthetic gasses mixing with the air or nitrous oxide.** Otherwise, it may cause risk.
- **If liquid splashes on the device or accessories, especially when liquids may enter the pipe or device, stop using and contact the service department.** Otherwise, it may cause risk.
- **Dispose of the packaging material in compliance with the applicable waste regulations and keep it away from children.** Otherwise, it may cause harm to the environment or children.
- **Please use approved accessories for the device and check that the device and accessories are working properly and safely before use.** Otherwise, the measurement result may be inaccurate or an accident may occur.
- **If the unit has accidentally become damp, it should be placed in a dry and ventilated place for some time to allow the moisture to dissipate.** Otherwise, the device may be damaged due to moisture.
- **Do not store and transport the device outside the specified environment.** Otherwise, it may cause measurement error.
- **It is recommended to check the device and accessories regularly for damage.** If you find any damage, discontinue use and contact the hospital biomedical engineer or our customer service immediately. **Do not disassemble, repair or modify the unit without authorization.** Otherwise, it cannot be accurately measured.
- **This device cannot be used on mobile transport platforms.** Otherwise, it may cause measurement error.
- **This device cannot be used on a tilted tabletop.** Otherwise, there is a risk of falling.
- **Dispose of packaging materials, waste batteries and end-of-life products in accordance with local laws and regulations. The end-of-life products and materials are properly disposed of by the user in accordance with the authority's decree.**
- Replacing accessories that were not provided by our company may lead to the occurrence of errors.
- **Without our company or other approved maintenance organizations trained service personnel should not try to maintain the product.**
- This device can only be used for one test object at a time.
- **If the small parts on the device are inhaled or swallowed, please consult**

a doctor promptly.

- The device and accessories are processed with allergenic materials. If you are allergic to it, stop using this product.
- **If the unit has a display error after pressing the power button, such as white screen, blurred screen or no display content, please contact our company.**
- The device shall comply with the standard IEC 80601-2-30: Particular requirements for basic safety and essential performance of automated non-invasive sphygmomanometers.

2.1 Operation for power adapter (not included)

Note

The device can be powered by a power adapter. Be sure to use the dedicated medical grade power adapter of this device. Otherwise, it may cause trouble

Dedicated power adapter must use AC 100 V~240 V. Otherwise, it may cause fire or electric shock.

When there is breakage of dedicated power adapter plug or wire, please do not use it. Otherwise, it may cause fire or electric shock.

Please do not plug or unplug the adapter on the socket with wet hands. Otherwise, it may cause electric shock or injury.

When using a power adapter to connect with the power socket, make sure the power socket is conveniently accessible, in order to timely disconnect from the power when emergency.

2.2 Operation for Battery

Note

Please use 4 "AA" size manganese or alkaline batteries, do not use batteries of other types. Otherwise, it may cause fire.

Do not mix old and new batteries and batteries of different types. Otherwise, it may cause battery leakage, heat, rupture, and damage to Electronic Sphygmomanometer.

Please do not confuse the positive and negative battery terminals. When the batteries power exhausts, replace with four new batteries at the same time.

Please take out the batteries when you do not use the device for a long time (3 months or more). Otherwise, it may cause battery leakage, heat, rupture, and damage to Electronic Sphygmomanometer.

If electrolyte of the batteries immodestly get in your eyes, immediately rinse with plenty of clean water.

It will cause blindness or other hazards, should immediately go to the nearest hospital for treatment.

If electrolyte of the batteries immodestly glues on the skin or the clothes, immediately rinse with plenty of clean water. Otherwise, it may hurt the skin.


Advice

- Do not strike or drop the device;
- Do not inflate before the cuff wraps around the arm;
- Do not inflect the cuff and the air tube forcibly.

The device applies measurement Blood Pressure (BP) and Pulse of adult.

Chapter 3 Main Unit

All products are in the box. Open the box and check if the product is complete.

- | | |
|------------------|----------------------|
| 1. Cuff Plug | 4. MEMORY Button |
| 2. Cuff Air Plug | 5. START/STOP Button |
| 3. Adult Cuff | 6. Display |

Adult Cuff + Extra-large adult cuff:

Specification: limb circumference 18-26 cm and 22-43 cm (middle part of upper arm), please choice suited cuff when measuring other.

Also included: User Manual

© The cuffs are consumables. Calculate by measuring 6 times a day (3 times each morning and evening), the service life of the cuff is about 1 year (using our experimental conditions).

© To measure the blood pressure correctly, please replace the cuff in time.

© If the cuff leaks, please contact our company to buy a new one. The cuff

purchased separately does not include the airway tube plug. When replacing, please do not throw the airway tube plug away, but install it on the new cuff.

Note

When the product and accessories described in this manual are about to exceed the period of use, they must be disposed according to relevant product handling specification. If you would like more information, please contact our company or representative organization.

Chapter 4 External Interfaces

Note

When removing NIBP cuff, please take the plug at the front of the windpipe to pull it out.

① Cuff socket ( is cuff identifier)

② Power adapter socket ( is power socket identifier)

Note

All analog and digital equipment connected to this device must be certified to IEC standards (such as IEC60950: Information technology equipment-Safety and IEC60601-1: Medical electrical equipment-Safety), and all equipment should be connected to in accordance with the requirement of the valid version of the IEC60601-1-1 system standard. The person connecting the additional equipment to the signal input and output port is responsible for whether the system complies with the IEC60601-1 standard.

Chapter 5 Battery/AC Adapter Installation

The product can be operated with batteries and AC adapter.

5.1 Battery Installation

- ① Demount the battery cover in the direction of the arrow.
- ② Install "AA" batteries according to polarities .

- ③ Slide to close the battery cover.

Icon “”: the batteries power will exhaust. Replace with four new batteries (of the same sort) at the same time. Test while low power may cause data deviation and other problems.

Turn the unit off before replacing the batteries.

Note

When the battery reaches the end of its life, or if the battery is found to have odor, deformation, discoloration or distortion, stop using the battery and dispose of the used battery in accordance with local regulations, otherwise it will cause environmental pollution.

5.2 Usage of power adapter

1. Connect the sphygmomanometer and the power adapter. Plug the power adapter plug into the power adapter socket on the back of the device
2. Please insert the power plug of the adapter into the AC 100 V~240 V socket.

Note

The device can be disconnected from the power supply network by unplugging the adapter plug.

When cut off the power supply, first cut off the connection of power socket and the regulated power supply, then cut off the connection of regulated power supply and the sphygmomanometer.

Please be sure to use dedicated medical grade power adapter.

Note

When regulated power supply and batteries are both used at the same time, the battery power will not be consumed.

Switch regulated power supply and battery as power supply when the device is off, otherwise, the device may shutdown due to power failure.

The device can be used normally after switching on, without waiting for the device to be prepared.

Chapter 6 Button Functions

6.1 Description for button operation

All the operations to the Electronic Sphygmomanometer are through buttons. The names of the buttons are written above them. They are:

- Left button is “M” button. Under “OFF” state, press this button to enter the review interface (refer to Chapter 8 for details).

- Right button is “START/STOP” button. Under “OFF” state, press this button to enter measurement mode, inflate the cuff to measure

blood pressure, press this button again to turn off the device.

6.2 Units setting

Under "OFF" state, press "M" button and "START/STOP" button simultaneously for 5 seconds to enter the setting interface, the default unit in this interface is "mmHg"; short press "M" button to switch the unit between "mmHg" and "kPa".

Note

■ The default unit of the device when leaving factory is mmHg.

Chapter 7 The Usage Method of Sphygmomanometer

7.1 Accurate Measurement Way

Measurement in quiet and relaxed state.

1. Adopt a comfortable sitting position, use back and arms to support the body.
2. Place your elbow on a table, the palm faces up and the body is relaxed.
3. The cuff is at the level of your heart.
4. Feet flat on the floor, and do not cross your legs.

Advice

Try to measure your blood pressure at the same time each day with the same arm and the same pose for constancy.

The high and low location of cuff will cause changes in measure results.

Do not touch the sphygmomanometer, cuff and windpipe during measure.

The measurements should be taken in a quiet place and the body should be relaxed.

Remain still 4~5 minutes before measurement.

Do not talk or move during the measurement. Relax the body, do not allow any muscle activity.

Wait 4~5 minutes between measurements.

Do not use precision instrument near the Sphygmomanometer.

Warning

During repeated measurements, the exact blood pressure value may not be measured due to congestion in the arm. Please measure when the blood flow is calm.

Repeated measurements over a longer period of time or limbs rubbing with the cuff may be accompanied by purpura, ischemia and nerve damage.

When measuring a patient, the color, warmth, and sensitivity of the distal limbs must be checked regularly. Once any abnormalities are observed, place the cuff in another position or stop the blood pressure measurement

immediately.

Please use the device at an environment of suitable temperature and humidity otherwise it will cause measurement errors.

Do not twist or wrap the airway tube. It can cause constant pressure in the cuff which can block blood flow and cause serious damage to the patient. Do not use the cuff on injured areas as this may cause more severe damage to the area.

Do not use the cuff in the area where the treatment is being performed inside blood vessel or the arteriovenous connection. This may cause temporary blockage of blood flow and cause injury to the patient.

Do not use the cuff on the side of the mastectomy.

When using the cuff to pressurize, some of the body's functions may temporarily weaken. Do not use the medical electrical measuring device in the corresponding arm position.

Do not move during measurement, it will have a delayed effect on the patient's blood flow.

The device needs to be placed for 2 hours from the minimum storage temperature to being ready for its intended use.

The device needs to be placed for 4 hours from the highest storage temperature to being ready for its intended use.

Note

The following conditions may also cause changes in the blood pressure measurement value.

If you take the measurement within one hour of eating, drinking alcohol or coffee, or smoking, exercising, bathing;

Using incorrect posture such as standing or lying down, etc;

The patient speaks or moves his body during measurement;

When measuring, the patient is nervous, excited, emotional instability;

The room temperature rises or falls sharply, or the environment of measurement often changes;

Measuring in a moving vehicle;

The high and low location of cuff will cause changes in measurement results;

Continuous measurement over a long period of time.

7.2 Applying the Cuff

Both left and right arm can be measured.

Bare your arm or wear tight-fitting clothing during the measurement.


Carry out the operation in a room with comfortable temperature.

When measuring, take off the thick clothing instead of rolling up the sleeves.

To measure as accurately as possible, make sure you put the cuff on

correctly (left arm).

- ① Insert the arm cuff air plug in the cuff socket of sphygmomanometer.
- ② Stretch the cuff into a barrel so that the arm can enter the barrel conformably.
- ③ The left arm is guided into the cuff, the air tube of the cuff passes the palm of the hand.
- ④ Place the cuff around your upper arm. Place the air tube on the inside of your forearm and align it using your middle finger.
- ⑤ The bottom of the cuff should be approximately 2cm~3cm above your elbow.
- ⑥ Once the cuff is tightly wrapped and fixed, there should be no gaps between the arm and the cuff.

7.3 Measurement BP

- ① Under "OFF" state, press "START/STOP" button to start measuring. Please keep the correct pose and still body state during the measurement, the body should not be moved. The "Movement" icon appears if patient moves and continuing the measuring may lead to inaccurate measurements.

If you want to abort the measurement

Press the 【START/STOP】 button, the device will stop inflating, and release the air from the cuff.

- ② The display of the measurement results appears after finishing the measurement. The pressure bar on the right side visually demonstrates the pressure level.

7.4 Confirm the measured value

① The World Health Organisation has established globally recognised standards for the assessment of hypertension (in the clinical setting).

Blood pressure level	Systolic pressure	Diastolic pressure
Normal	pressure < 130mmHg	pressure < 85mmHg
Normal systolic value	130mmHg ≤ pressure ≤ 139mmHg	85mmHg ≤ pressure ≤ 89mmHg
Mild hypertension	140mmHg ≤ pressure ≤ 159mmHg	90mmHg ≤ pressure ≤ 99mmHg
Moderate hypertension	160mmHg ≤ pressure ≤ 179mmHg	100mmHg ≤ pressure ≤ 109mmHg
Severe hypertension / high blood pressure	180mmHg ≤ pressure	110mmHg ≤ pressure

② Pressure bar on the right side

A SYS value of more than 135 mmHg or a DIA value of more than 85 mmHg is used as criteria for high blood pressure (in the home environment) and the pressure bar on the right-hand side lights up red.

A SYS value below 135 mmHg and a DIA value below 85 mmHg are used as criteria for normal blood pressure and the pressure bar on the right side lights up solid green.

The number of lights represents the blood pressure range.

*Diagnosing and treating yourself based on measurement results can be dangerous. Follow the instructions of your doctor.

Note

Wait at least 4-5 minutes between measurements.

- With repeated measurements, the exact blood pressure value may not be measured due to congestion in the arm. Please measure when the blood flow is calm.
- When the screen displays “Err”, the measure can't be carried out correctly.
- Irregular pulse icon is displayed in the measurement results if the pulse interval is irregular during measuring, which may cause it is unable to take measurement correctly. Please keep quiet and remeasure. If the irregular pulse icon appears frequently, please consult a doctor.
- The minimum value of the patient's physiological signal is the minimum limit that the device can measure. The device may obtain inaccurate measurement results when it is operated below the minimum amplitude or minimum value of the patient's physiological signal.

*The device switches off automatically after five minutes without operation, even if you have forgotten to switch off the unit.

Chapter 8 Memory function

The device can store NIBP values automatically and display up to 199 set of measurement results.

If 199 measurement data sets have been stored in device, the oldest data set is overwritten when the 200th data set is saved. If there are no measurement values, the memory values cannot be numbered.

Memory function cannot be used during measuring.

When there are no measurement values, “---” will display on the review interface.

8.1 Review the Memory Value

Under “OFF” state, press the “M” button to display the average value of the latest three set of data. If the number of measurement data is less than three groups, it is automatically supplemented. Continue to press “M” button in current interface to view all measurement records.

8.2 Delete Memory Values

1. Users can delete all memory values of the current user instead of separately delete one memory value.

2. Press the “M” key and the “START/STOP” key simultaneously for more than 5 seconds under the memory interface, as soon as “DEL” appears on the screen, all memory values are deleted.

Caution

When querying the measurement records, please press “M” button continuously to query one by one.

Chapter 9 Symbols

Your device may not contain all the following symbols.

Signal	Description	Signal	Description
	Attention! Please refer to the accompanying document (the user manual).		Attention! Please refer to the accompanying document (the user manual).
SYS	Systolic pressure	DIA	Diastolic pressure
MAP	Mean blood pressure	PUL	Pulse rate (bpm)
IPX0	Enclosure protection grade	EMC	Electromagnetic compatibility
	Recyclable	P/N	Material code of manufacturer
	Batch code		Use by date
	This way up		Fragile, handle with care

	Keep dry		Storage atmospheric pressure limitation
	Storage temperature limitation		Storage humidity limitation
	Manufacturer		Date of manufacture
	Batteries Power	SN	Serial number
	Inflating		Deflating
	Waste disposal mark, this symbol indicates that the waste of electrical and electronic equipment cannot be disposed as an unclassified municipal waste and must be recovered separately.		This item is compliant with Medical Device Directive 93/42/EEC of June 14, 1993, a directive of the European Economic Community.
	Class II equipment		Type BF applied parts
	European Representative		Irregular pulse
	Socket for power adapter		Interface for connecting cuff
	Voice closed		Voice enabled
	Large movement during measurement		Cuff tied properly
	Artery indicator label		Swiss authorized representative
	UK Authorized Representative		UKCA marking

Chapter 10 Error Message

If "Err" appears at the high pressure position and the error number appears at the low pressure position, the measurement is not normal. (Error numbers are 02,04,06-16,19)

Error Mark	Causes	Solutions
Err02 Err15	Function abnormal	Please contact us
Err04	Low battery	Please replace the battery or link the adapter
Err06	The cuff is not wrapped	Wrap the cuff correctly (refer to

	correctly.	Chapter 10)
Err07	Cuff plug fall off	Make sure the cuff plug is securely inserted in the windpipe (refer to Chapter 10)
Err08	Air pressure error	Keep arm and body still, measure again
Err09	The pulse signal is too weak or the cuff is loose.	Wrap the cuff correctly (refer to Chapter 10)
Err10	Cuff is blocked or squeezed	Wrap the cuff correctly (refer to Chapter 10)
Err14	Cuff leakage	Replace with a new cuff
Err11 Err12 Err13	The signal amplitude is too big owing to the arm or body moving or other reasons when measuring	Keep arm and body still, measure again
Err16 Err19	It takes too much time	

Chapter 11 Troubleshooting

Abnormal Phenomenons	Causes	Solutions
BP measurement values too high or too low.	Cuff is not connected correctly.	Correctly connect cuff.
	Spoken or arm moved during measurement	Keep quiet and restart a measurement.
	The arm is squeezed too tightly.	Take off the clothes, and restart a measurement
No pressure	Cuff leakage	Buy a new cuff.
	The cuff windpipe is not correctly connected with cuff	Correctly connect.
	Cuff does not inflate	Contact us.
Cuff deflates in a short time	Loose cuff	Correctly tangle cuff.
The measurement cannot be continued even if the measurement button is pressed.		Switch the unit on again and start a new measurement.
Abrupt switching off in adding pressure	No use for a long time, the batteries can be exhausted owing to the changed temperature	Replace all four batteries with new ones.

On/Off button is held down but the unit does not start up	Batteries are worn	Replace all four batteries with new ones.
	The battery polarities are reversed	Check the battery installation for proper placement of the battery polarities.
Cuff inflation starts before pressing the measurement button		Stop using the device and contact us.
Cuff never deflates		Stop using the device and contact us.
Air pressure error	Deflation error	Pull out the cuff to deflate. Stop using the device and contact us.
	Others	Keep arm and body still, measure again.
No pressure value is displayed or the value remains unaltered when cuff inflating		Pull out the cuff to deflate. Stop using the device and contact us.
Other phenomenon		Switch on the power once again and restart an operation. Replace the batteries. If this does not help, please contact us.

Chapter 12 Maintenance, Cleaning and Keeping

*Please do obey the precautions and correct operating methods in this user manual. Otherwise, we will not be responsible for any fault.

Warning

Remove the batteries before cleaning. The accessories and the main unit must be separated for cleaning.

Maintenance is not allowed during device using.

Do not squeeze the rubber tube on the cuff.

Caution

- High pressure disinfection to the device and accessories is not allowed.
- Do not let water or cleaning agent flow into the socket to avoid device damage.
- Do not soak the device and accessories in liquid.
- If any damage or deterioration of the device and accessories is found, please do not use it.

Maintenance:

- Clean the device and accessories regularly. It is recommended to clean

them once a month. When the device or accessory gets dirty, use a dry and soft cloth to wipe. If they are very dirty, it is available to dip the soft cloth into water or mild detergent and wring it out before using it for cleaning.

- The device shall be inspected and calibrated regularly (or according to inspection standard of hospital). The inspection can be carried out in appointed institutions, or by professional personnel or contact us for inspection. Under the setting interface, Press the “START/STOP” button once, after “CAL” appears on the screen, press “M” button for more than 15 seconds to enter the static pressure interface.

 Advice

- Do not use gasoline, volatile oil, diluent, etc. to wipe the device.
- Do not clean or wet the cuff.

Storage:

 Advice

- Do not expose the device in direct sunlight for long time, otherwise the display screen maybe damaged.
- The basic performance and safety of the device are not affected by the dust or cotton wool in home environment, while the device shall not be placed in locations with high temperature, humidity or dust.
- Aged cuff may result in inaccurate measurement; please replace the cuff periodically according to the user manual.
- To avoid device damage, keep the device out the reach of children and pets.
- Avoid placing the appliance near extremely high temperatures, such as a fireplace, as this may affect the performance of the device.
- Do not store the device with chemical medicine or corrosive gas.
- Do not place the unit near water.
- Do not place the unit in a location where it is subject to inclination, vibration or shock.
- Take the batteries out if the device is not to be used for three months or longer.

Chapter 13 NIBP Specification

Name	Electronic Sphygmomanometer
The degree of protection against ingress of water	IPX1
Display	LED
NIBP Specifications	
Measurement Method	Oscillometric method

Working mode		Automatic
Operation mode		Continuous operation
Pressure Range		0~297 mmHg (0~39.6 kPa)
Measurement range	Pressure	SYS : 30~270 mmHg (4~36 kPa) DIA : 10~220 mmHg (1.3~29.3 kPa)
	Pulse: 40~240/min	
Inflation		160±5 mmHg (21.33±0.67 kPa)
Overpressure protect		297±3 mmHg (39.6±0.4 kPa)
Resolution		Pressure: 1 mmHg (0.133 kPa)
Accuracy		Static pressure: ±3 mmHg (±0.4 kPa)
Error	The BP value measured by the device is equivalent with the measurement value of Stethoscopy, perform clinical verification in accordance with the requirements in ISO 81060-2: 2013, whose error meets the followings: Maximum mean error: ±5 mmHg Maximum Standard deviation: 8 mmHg	
Operating Temperature/ Humidity		+5°C~40 °C. 15%RH~85%RH (no condensation)
Transport	Transport by general vehicle or according to the order contract, avoid pounded, shake and splash by rain and snow in transportation.	
Storage	Temperature: -20 °C~+55 °C; Relative humidity: ≤95 %; No corrosive gas and drafty.	
Atmospheric pressure	700 hPa~1060 hPa	
Power supply	4 "AA" alkaline batteries, AC Adapter (AC, 100 V-240 V, optional)	
Rated current	≤ 600 mA	
Battery life	When the temperature is 23 °C, limb circumference is 270 mm, the measured blood pressure is normal, 4 "AA" alkaline batteries cab be used about 300 times.	
Main Unit Dimensions	129*101*72 mm	
Main Unit Weight	281 gram (without batteries)	
Safety classification	Class II equipment (power supplied by power adapter)/Internally powered equipment (power supplied by batteries), Type BF applied part	

Service life	The service life of the device is five years or 10000 times of BP measurement.
Date of manufacturer	See label
Accessories	Standard Configure: Adult Cuff: limb circumference 22-32 cm (upper arm center), User Manual, four "AA" alkaline batteries AC Adapter: Input: voltage: AC 100 V~240 V frequency: 50 Hz/60 H Rated current: AC 150 mA Output: DC 5.0 V $\pm$ 0.2 V 1.0 A

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CONTEC MEDICAL SYSTEMS CO, LTD

ADD: No.112 Qinhuang West Street, Economic & Technical Development Zone, Qinhuangdao, Hebei Province
066004, P.R. China

Prolinx GmbH

ADD: Brehmstr. 56, 40239 Duesseldorf, Germany
Tel: 0049 211 3105 4698
E-mail: med@eulinx.eu

MedPath Limited

ADD: 27 Old Gloucester Street, London WC1N 3AX, United Kingdom,
E-mail: UK@medpath.pro, Web: www.medpath.pro
Tel: +44 776 280 3783

Tinovamed GmbH

ADD: Zielmattenring 6, 4563 Gerlafingen, Switzerland
info@tinovamed.ch, tinovamed.ch
+41 32 566 71 72

Bei Fragen wenden Sie sich bitte an/If you have any questions, please contact:/Si tiene alguna pregunta, póngase en contacto con nosotros/Si vous avez une question s'il vous plaît contacter/Se si dispone di una domanda si prega di contattare/ Als u vragen heeft, kunt u contact opnemen met/ W razie jakichkolwiek pytań prosimy o kontakt/ Om du har några frågor, vänligen kontakta:

Novidion GmbH, Fuggerstr. 30, 51149 Köln, Germany

Tel.: +49 2203 - 9885 200, Fax: +49 2203 - 9885 206

info@pulox.de

www.pulox.de

WEEE: DE24355330

Sollten Sie beim Lesen dieser Anleitung Probleme haben, so können Sie diese auch auf unserer Website downloaden/ In case you have any problems reading this manual, you can download it from our website/ En caso de que tenga algún problema para leer este manual, puede descargarlo de nuestro sitio web/ Si vous rencontrez des problèmes pour lire ce manuel, vous pouvez le télécharger depuis notre site Web/ In caso di problemi nella lettura di questo manuale, è possibile scaricarlo dal nostro sito Web/ Als u problemen hebt met het lezen van deze instructies, kunt u ze ook downloaden van onze website/ Jeśli masz problemy z odczytaniem tej instrukcji, możesz ją również pobrać z naszej strony internetowej/ Om du har problem med att läsa dessa instruktioner kan du också ladda ner dem från vår webbplats:

<https://www.pulox.de/Anleitungen-Datenblaetter-Retouren>

Versionsnummer: 2.04 // 30.10.2023

Appendix

Guidance and manufacturer's declaration – electromagnetic emissions-for all EQUIPMENT and SYSTEMS

Guidance and manufacturer's declaration – electromagnetic emission		
The device is intended for use in the electromagnetic environment specified below. The customer of the user of the device should assure that it is used in such and environment.		
Emission test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emission CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC61000-3-2	Class A	
Voltage fluctuations/ flicker emissions	Complies	

IEC61000-3-3		
--------------	--	--

**Guidance and manufacturer's declaration – electromagnetic immunity –for
all EQUIPMENT and SYSTEMS**

Guidance and manufacturer's declaration – electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floor are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines	±2kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T)for 0.5 cycle 40% U_T (60% dip in U_T)for 5 cycles 70% U_T (30% dip in U_T)for 25 cycles <5% U_T (>95% dip in U_T)for 5 sec	<5% U_T (>95% dip in U_T)for 0.5 cycle 40% U_T (60% dip in U_T)for 5 cycles 70% U_T (30% dip in U_T)for 25 cycles <5% U_T (>95% dip in U_T)for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. The device can continue the operation during power mains interruptions due to the usage of battery.

Power frequency (50/60Hz) magnetic field IEC61000-4-8	3A/m	3A/m	Mains power quality should be that of a typical commercial or hospital environment.
---	------	------	---

NOTE: U_T is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration – electromagnetic immunity –for EQUIPMENT and SYSTEMS that are not LIFE-SUPPORTING

Guidance and manufacturer's declaration – electromagnetic immunity

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
			Portable and mobile RF communications equipment should be used no closer to any part of the <i>device</i> , including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Conducted RF IEC6100 0-4-6	3 V/m 150 kHz to 80 MHz	3 V/m	Recommended separation distance $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$ $d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$ 80 MHz to 800 MHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$ 800 MHz to 2.5 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b
Radiated RF IEC6100 0-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	

			Interference may occur in the vicinity of equipment marked with the following symbol: 
--	--	--	---

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

A Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.

B Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the EQUIPMENT or SYSTEM – for EQUIPMENT or SYSTEM that are not LIFE-SUPPORTING

Recommended separation distances between portable and mobile RF communications equipment and the <i>device</i>			
The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter(m)		
	150 kHz to 80 MHz $d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74

1	1.17	1.17	2.33
10	3.69	3.69	7.38
100	11.67	11.67	23.33

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Warning

- Active medical devices are subject to special EMC precautions and they must be installed and used in accordance with these guidelines.
- Electromagnetic fields can affect the performance of the device, so other equipment used near the equipment must meet the appropriate EMC requirements. Mobile phones, X-rays, or MRI devices are possible interference sources, as they emit high-intensity electromagnetic radiation.
- The use of ACCESSORIES, transducers and cables other than those specified, with the exception of transducers and cables sold by the MANUFACTURER of the the device as replacement parts for internal components, may result in increased EMISSIONS or decreased IMMUNITY of the ME EQUIPMENT or ME SYSTEM.
- The device should not be used when they are close to or stacked with other equipment, if necessary, please observe and verify that they can operate normally in the configurations.
- Devices or systems may still be interfered by other equipment, even if other equipment meets the requirements of the corresponding national standard.
- The device requires special precautions for electromagnetic compatibility (EMC) and requires qualified personnel to install and use in accordance with the EMC information provided below.
- The device should not contact the pins of connectors marked with an ESD warning symbol, unless electrostatic discharge precautions are used, the device should not connect to these connectors.
- In order to avoid the accumulation of electrostatic charge, it is

recommended to store, maintain and use the equipment at a relative humidity of 30 % or more. The floor should be covered with ESD dissipated carpets or similar materials. In the use of the components, non-synthetic clothing should be worn.

- In order to prevent electrostatic discharging to the ESD-sensitive parts of the device, the personnel should contact the metal frame of the components or the large metal objects near the device. When using the device, especially when it is possible to contact the ESD-sensitive parts of the device, the operator should wear a grounded bracelet designed for ESD-sensitive devices. For more information on proper use, please refer to the instructions provided with the bracelet.
- All potential users are advised to understand the ESD warning symbols and receive training on ESD precautions.
- The most basic content of the ESD precautionary procedure training should include an introduction to electrostatic charge physics, voltage level in the conventional case, and damage to the electronic components when the operator with electrostatic charge contacts them. In addition, the methods for preventing electrostatic buildup, and the manner and reasons for the release of human body static electricity to the ground or equipment frame or the use of a bracelet to connect the human body to the equipment or the ground before establishing the connection should be described.

The following cable types must be used to ensure that they comply with interference radiation and immunity standards:

Name	Length (m)
Power adapter cable	1.5