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Plastic cupping glasses set with vacuum pump 19 pcs.

Instruction manual

Dear Customer,

we are pleased that you have chosen a product from our range.

Please read these operating instructions carefully and completely before use and keep them accessible for other users. If you have any questions, always consult your attending physician before using the product. Follow the instructions.

With kind regards

Your Novidion Team

Mode of action

Cupping is a traditional therapy method that has been used since ancient times. By sucking the skin and the resulting vacuum in the cupping glass, a negative pressure is created which stimulates the acupuncture points on the surface of the body. Cupping can be used to stimulate blood circulation and metabolism, to drain blood stasis, to prevent cellulite and to tighten connective tissue.

Safety instructions

1. The following parts of the body must not be cupped: Surroundings of the heart, on large blood vessels, eyes, ears, nose, nipples and in case of skin diseases or acute inflammation of the affected skin area.
2. Do not apply the cupping glasses on open, swollen, bruised or purulent skin areas.
3. Contraindications in diseases with bleeding tendency, renal or cardiac insufficiency, acute bone fracture, anemia, critical hypertension, neurodermatitis.
4. Cupping glasses should not be used by pregnant women.
5. To prevent the cupping glasses from falling down, the patient should not move during the treatment.
6. If the suction marks do not disappear after the first cupping, the cupping glasses must not be applied in the same place again.
7. If the skin turns red or itches after removing the cupping glasses, apply a special oil or cream to reduce rashes after cupping.
8. During the treatment, observe the patient continuously and remove the cupping glasses in case of pain, discomfort or unusual muscle movements.
9. Always treat the product with care and hygiene, make sure that you clean and disinfect it thoroughly after each application (see chapter: Cleaning). This is the only way to avoid cross-infections.
10. Do not drop the product, it could break. Do not store or use it in an environment with high humidity or high temperature.
11. Do not use the device in an explosive or hazardous environment.

12. Keep the cupping glasses away from children. The child may get hurt or damage the cupping glass.

Application

1. Choose the place you want to cupping, clean the skin place and remove the hair in case of hairy skin. Apply some grease, such as massage creams or lotions, to the skin area before use.
2. Select the size of the cupping glass according to the skin area. The larger the skin site is, the larger the cupping glass should be.
3. Loosen the air valve at the top of the cupping and make sure that air can flow through it.
4. Place the cupping glass on the corresponding skin area so that the rim of the cupping glass is tight and without gaps on the skin.
5. Place the vacuum pump on top of the air valve. Pull the handle of the vacuum pump several times so that a vacuum is created and the skin is slightly sucked by the cupping glass. The vacuum in the cupping glass can be regulated by pumping. Normally, the skin in the cupping glass is lifted about 5-15 mm.
6. Hold the cupping glass with one hand and remove the vacuum pump from the cupping glass. Press the air valve on the top of the cupping glass to prevent air from escaping. A cupping application usually takes about 10-15 minutes. However, the treatment duration should be adjusted to the individual patient and shortened if necessary.
7. To remove the cupping glass, hold it firmly with one hand. With the other hand, gently pull up the air valve so that air flows into the cupping glass and the negative pressure decreases until it disappears. Now you can detach the cupping glass from the skin area.

Symbol	Description
	See operating instructions
SN	Serial number
MD	Medical device
	Manufacturer
EC REP	European representative
	Importer
CE	Product complies with Directive 2014/30/EU of the European Parliament and of the Council of February 26, 2014 and Directive 2001/95/EC of the European Parliament and of the Council of December 3, 2001.

Application with the cupping hose

If you want to cupping yourself on hard-to-reach parts of the body, such as the back, neck, back of the legs, etc., you can use the cupping hose. Proceed in the same way as in the point "Application", except that you connect the cupping glass and the vacuum pump via the cupping hose.

Maintenance and upkeep

Always treat the product with care and hygienically. Clean it thoroughly after each use and disinfect it with medical alcohol. This is the only way to avoid cross-infections.

Do not use damaged cupping glasses. Check the cupping glasses before each use.

Store the cupping glass set away from dust and moisture.

Warranty and service

We provide 1 year warranty for material and manufacturing defects of the product. The warranty does not apply:

- in case of damage caused by improper operation
- for wearing parts
- for defects that were already known to the customer at the time of purchase
- in the case of the customer's own fault

The legal warranties of the customer remain unaffected by the warranty. For the assertion of a warranty claim within the warranty period, the customer must provide proof of purchase. The warranty must be claimed within a period of 1 year from the date of purchase from:

Novidion GmbH, Fuggerstr. 30, 51149 Cologne, Germany

Further rights are not granted to the customer due to the warranty.

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Novidion GmbH, Fuggerstr. 30, 51149 Köln, Tyskland
Jiangmen Xinli Medical Apparatus and Instruments Co.,Ltd
Lotus NL B.V.