HC-10

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**Operating Instructions (EN)**

Dear customer,

we are pleased that you have decided to purchase a product from our range. Please read these operating instructions carefully and keep them accessible for other users. Please pay attention to the instructions.

With kind regards

Your Novidion Team

Instructions and warnings for gel gloves**Cold Therapy:**

1. Place the gloves flat in the freezer for 90 minutes.
2. Take the gloves out of the freezer and test the temperature carefully according to your personal feeling.
3. If the temperature is acceptable, you can apply the gloves.
4. If the temperature is too cold, remove the gloves and wait until the temperature is comfortable.

Heat Therapy:

1. Spread the gel evenly in the gloves and place the gloves flat in the center of the clean microwave turntable (without touching the sides). Make sure that all additional microwave functions (such as "browning") are disabled.
2. At maximum wattage, heat the microwave at full power for 30 seconds. Then turn the gloves over and heat for another 30 seconds.
3. Remove the gloves from the microwave. Knead well to distribute heat evenly and gently check to see if desired temperature is reached. (The temperature will continue to rise slightly).
4. If the gloves feel too hot, remove them and put them back on when the temperature is more comfortable.
5. For additional warmth, place back in microwave and heat at 30 sec intervals.
6. Watch the gloves as they heat, as extreme heating can cause them to tear or leak.
7. When the temperature is comfortable, put on the gloves and remain seated during use.

Note: Do not apply cold or heat therapy more than 20 minutes a time.

Instructions for Care: Gently blot gloves with warm water and air dry.

ATTENTION:

These gloves are only suitable for cold & heat therapy.

CAUTIONS:

- Individuals with circulatory or diabetic issues should consult a physician before use
- Do not apply excessive weight to gloves by standing or sitting on them
- External use only
- Unattended use of gloves by children, elderly or incapacitated persons is dangerous
- Monitor for insensitivities to product temperature
- Do not apply to open wounds
- Discard if punctured or leaking
- If contents are swallowed, drink several glasses of water and contact Poison Control Center immediately
- If contents are exposed to skin or eyes, rinse thoroughly and consult a physician
- Do not use product with any ointments, creams or topical pain relievers
- If condition persists, consult a physician
- Keep out of reach of children

Disposal: Please dispose of the packaging according to local regulations.

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Novidion GmbH
Wuhan Pakcare Bio-technology Co., LTD.
Share Info GmbH